

Product datasheet for **TP314343**

Aminomethyltransferase (AMT) (NM_000481) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human aminomethyltransferase (AMT), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC214343 representing NM_000481 Red =Cloning site Green =Tags(s)
	<p>MQRAVSVVARLGFRLQAFPPALCRPLSCAQEVLRRTPLYDFHLAHGGKMVAFAGWVSLPVQYRDSHTDSDL HTRQHCSLFDVSHMLQTKILGSDRVKLMESLVGDIARPNQGTLSLFTNEAGGILDDLIVTNTSEGLH YVSNAGCWEKDLALMQDKVRELQNQGRDVGLEVLNALLALQGPTAAQVLQAGVADDLRKLPMFMTSAVM EVFGVSGCRVTRCGYTGEDGVEISVPVAGAVHLATAILKNPEVKLAGLAARDSLRLEAGLCLYGNDIDEH TTPVEGSLSWTLGKRRRAAMDFFGAKVIVPQLKGRVQRRRVGLMCEGAPMRAHSPILNMEGTKIGTVTSG CPSPSLKKNVAMGYVPCYSRPGTMLLVEVRRKQQMAVVSKMPFVPTNYYTLK</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	43.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_000472
Locus ID:	275



[View online »](#)

UniProt ID:	<u>P48728, A0A024R2U7</u>
RefSeq Size:	2117
Cytogenetics:	3p21.31
RefSeq ORF:	1209
Synonyms:	GCE; GCST; GCVT; NKH
Summary:	This gene encodes one of four critical components of the glycine cleavage system. Mutations in this gene have been associated with glycine encephalopathy. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2011]
Protein Pathways:	Glycine, serine and threonine metabolism, Metabolic pathways, Nitrogen metabolism, One carbon pool by folate

Product images:



Coomassie blue staining of purified AMT protein (Cat# TP314343). The protein was produced from HEK293T cells transfected with AMT cDNA clone (Cat# [RC214343]) using MegaTran 2.0 (Cat# [TT210002]).