

## Product datasheet for **TP314360M**

### PCDHGA6 (NM\_018919) Human Recombinant Protein

#### Product data:

**Product Type:** Recombinant Proteins  
**Description:** Purified recombinant protein of Homo sapiens protocadherin gamma subfamily A, 6 (PCDHGA6), transcript variant 1, 100 µg  
**Species:** Human  
**Expression Host:** HEK293T  
**Expression cDNA Clone or AA Sequence:** >RC214360 representing NM\_018919  
**Red**=Cloning site **Green**=Tags(s)

MAPPQRHPQRSEQVLLLTLLGTLWGAAAAQIRYSIPEELEKGSFVGNIVKDLGLEPQELAEHGVRIVSRG  
RMQLFSLNPRNGSLVTAGRIDREELCAQSPRCLVSFNILVEDKLNLYPVEVEIVDINDNTPRFLKEELEV  
KILENAAPSSRFPLMEVYDPDVGMNSLQGFKLSGNSHFSVDVQSEAHGPKYPELVLEGLDREGEAVYRL  
VLTAMDGGDPVRSVAQILVTVLDVNDNTPMFTQPVYRVSVPENLPVGTTPVLAVTATDQDEGVHGEVTYS  
FVKITEKISQIFCLNVLGTGEISTSANLDYEDSSFYELGVEARDGPGLRDRAKVLITLDVNDNVPEVVT  
SGSRTIAESAPPGTVIALFQVFDNRDSSLNGLVTCISIPRSLPFELEKSVGNYYRLVTNAALDREEVFLYNI  
TVTATDKGTPLSTETIISLNVADTNDNPPTFPHSSYSVYVLENNPRGASIFSVNALDPDQNAQVSY  
LAEDTLQGAPLSSYVSINSDTGILYALRSFDYEQRLDLQWWTASDSGDPPLSSNVLSLFLVDQNDNAP  
EILYPALPTDGGSTGVELAPRSAEPGYLVTKVAVDRDQSGQNAWLSYRLLKASEPGLFSVGLHTGEVTRAR  
ALLDRDALKQSLVAVQDHGQPPLSATVTLTAVADRIDLADLGSLEPSAKPNDSDLTLYLVVAVAAV  
SCVFLAFVIVLLALRLQRWHKSRLQASGGGLASMPGSHFVGVGVRFLQYTSHEVSLTADSRKSHLIF  
PQPNYADTLINQESYEKSEPLLITQDLLETGKGEPRQLQAPPNTDWRFSQAQRPGTSGSQNGDDTGTWPN  
NQFDTEMLQAMILASASEAADGSSTLGGGAGTMGLSARYGPQFTLQHVPDYRQNVYIPGSNATLTNAAGK  
RDGKAPAGGNGNKKKSGKKEKK


**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

**Tag:** C-Myc/DDK  
**Predicted MW:** 97.7 kDa  
**Concentration:** >0.05 µg/µL as determined by microplate BCA method  
**Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining  
**Buffer:** 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol



[View online »](#)

<b>Preparation:</b>	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
<b>Note:</b>	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
<b>Storage:</b>	Store at -80°C.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_061742</a>
<b>Locus ID:</b>	56109
<b>UniProt ID:</b>	<a href="#">Q9Y5G7</a>
<b>RefSeq Size:</b>	4605
<b>Cytogenetics:</b>	5q31.3
<b>RefSeq ORF:</b>	2796
<b>Synonyms:</b>	PCDH-GAMMA-A6
<b>Summary:</b>	<p>This gene is a member of the protocadherin gamma gene cluster, one of three related clusters tandemly linked on chromosome five. These gene clusters have an immunoglobulin-like organization, suggesting that a novel mechanism may be involved in their regulation and expression. The gamma gene cluster includes 22 genes divided into 3 subfamilies. Subfamily A contains 12 genes, subfamily B contains 7 genes and 2 pseudogenes, and the more distantly related subfamily C contains 3 genes. The tandem array of 22 large, variable region exons are followed by a constant region, containing 3 exons shared by all genes in the cluster. Each variable region exon encodes the extracellular region, which includes 6 cadherin ectodomains and a transmembrane region. The constant region exons encode the common cytoplasmic region. These neural cadherin-like cell adhesion proteins most likely play a critical role in the establishment and function of specific cell-cell connections in the brain. Alternative splicing has been described for the gamma cluster genes. [provided by RefSeq, Jul 2008]</p>
<b>Protein Families:</b>	Transmembrane

**Product images:**

Coomassie blue staining of purified PCDHGA6 protein (Cat# [TP314360]). The protein was produced from HEK293T cells transfected with PCDHGA6 cDNA clone (Cat# [RC214360]) using MegaTran 2.0 (Cat# [TT210002]).