

Product datasheet for **TP314408L**

Integrin alpha 4 (ITGA4) (NM_000885) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human integrin, alpha 4 (antigen CD49D, alpha 4 subunit of VLA-4 receptor) (ITGA4), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC214408 representing NM_000885 Red =Cloning site Green =Tags(s)

MAWEARREPGPRRAAVRETVMLLLCLGVPTGRPYNVDTESALLYQGPHTLFGYSVLHSHGANRWLLVG
 APTANWLANASVINPGAIYRCRIGKNPGQTCEQLQLGSPNGEPCGKTCLEERDNQWLGVTLSRQPGENG
 IVTCGHRWKNIFYIKNENKLP TG GCGYVPPDLRTELSKRIAPCYQDYVKKFGENFASCQAGISSFYTKDL
 IVMGAPGSSYWTGSLFVYNITTNKYKAFLDKQNQVKFGSYLGYSVGAGHFRSQHTTEVWGAPQHEQIGK
 AYIFSIDEKELNILHEMKGKGLGSYFGASVCAVDLNADGFSDLLVGAPMQSTIREEGRV FVYINSGSGAV
 MNAMETNLVGS DKY AARFGESIVNLGDIDNDGFEDVAIGAPQEDDLQGAIIYNGRADGISSTFSQRIEG
 LQISKLSLMFGQSISGQIDADNNGYVDVAVGAFRSDSAVLLRTRPVVIDASLSHPESVNRTKFDCVEN
 WPSVCIDLTLCSYKGEVPGYIVLFYNMSLDVNRKAESPFRFYSSNGTSDVITGSIQVSSREANCRTH
 QAFMRKDVRDILTPIQIEAAYHLGPHVISKRSTEEFPLQPILOKKEKDIMKKTINFARFAHENC SAD
 LQVSAKIGFLKPHENKTYLAVGSMKTLMLNVSLFNAGDDAYETTLHVKLPVGLYFIKILEEEKQINCEV
 TDNSGVVQLDCSIGIYVDHLSRIDISFLLDVSSLSRAEEDLSITVHATCENEEEMDNLKHRSRVTAIPL
 KYEVKLT VHGFVNPTSFVYGSNDENEPETCMVEKMNLTFHVINTGNSMAPNVSVEIMVPNSFSPQTDKLF
 NILDVQTTTGECHFENYQRVCALEQQKSAMQTLKGIVQFLSKTDKRLLYCIKADPHCLNFLCNFGKMESG
 KEASVHIQLEGRPSILEMDETSALKFEIRATGFPEPNRVIELNKDENVAHVLLLEGLHHQRPKRYFTIVI
 ISSLLLGLIVLLLISYVMWKAGFFKRQYKXILQEENRRDSWSYINSKSNDD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	111.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol



[View online »](#)

Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_000876
Locus ID:	3676
UniProt ID:	P13612
RefSeq Size:	6082
Cytogenetics:	2q31.3
RefSeq ORF:	3096
Synonyms:	CD49D; IA4
Summary:	The gene encodes a member of the integrin alpha chain family of proteins. Integrins are heterodimeric integral membrane proteins composed of an alpha chain and a beta chain that function in cell surface adhesion and signaling. The encoded preproprotein is proteolytically processed to generate light and heavy chains that comprise the alpha 4 subunit. This subunit associates with a beta 1 or beta 7 subunit to form an integrin that may play a role in cell motility and migration. This integrin is a therapeutic target for the treatment of multiple sclerosis, Crohn's disease and inflammatory bowel disease. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2015]
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Transmembrane
Protein Pathways:	Arrhythmogenic right ventricular cardiomyopathy (ARVC), Cell adhesion molecules (CAMs), Dilated cardiomyopathy, ECM-receptor interaction, Focal adhesion, Hematopoietic cell lineage, Hypertrophic cardiomyopathy (HCM), Leukocyte transendothelial migration, Regulation of actin cytoskeleton

Product images:



Coomassie blue staining of purified ITGA4 protein (Cat# [TP314408]). The protein was produced from HEK293T cells transfected with ITGA4 cDNA clone (Cat# [RC214408]) using MegaTran 2.0 (Cat# [TT210002]).