

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TP314598

NTH1 (NTHL1) (NM_002528) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human nth endonuclease III-like 1 (E. coli) (NTHL1), 20 μg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC214598 representing NM_002528 <mark>Red</mark> =Cloning site Green=Tags(s)
	MCSPQESGMTALSARMLTRSRSLGPGAGPRGCREEPGPLRRREAAAEARKSHSPVKRPRKAQRLRVAYEG SDSEKGEGAEPLKVPVWEPQDWQQQLVNIRAMRNKKDAPVDHLGTEHCYDSSAPPKVRRYQVLLSLMLSS QTKDQVTAGAMQRLRARGLTVDSILQTDDATLGKLIYPVGFWRSKVKYIKQTSAILQQHYGGDIPASVAE LVALPGVGPKMAHLAMAVAWGTVSGIAVDTHVHRIANRLRWTKKATKSPEETRAALEEWLPRELWHEING LLVGFGQQTCLPVHPRCHACLNQALCPAAQGL
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	34.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP 002519</u>
Locus ID:	4913




View online »

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	NTH1 (NTHL1) (NM_002528) Human Recombinant Protein – TP314598
UniProt ID:	<u>P78549, E5KTI5</u>
RefSeq Size:	1097
Cytogenetics:	16p13.3
RefSeq ORF:	936
Synonyms:	FAP3; hNTH1; NTH1; OCTS3
Summary:	The protein encoded by this gene is a DNA N-glycosylase of the endonuclease III family. Like a similar protein in E. coli, the encoded protein has DNA glycosylase activity on DNA substrates containing oxidized pyrimidine residues and has apurinic/apyrimidinic lyase activity. [provided by RefSeq, Oct 2008]
Protein Families:	Druggable Genome, Stem cell - Pluripotency
Protein Pathway	s: Base excision repair

Product images:



Coomassie blue staining of purified NTHL1 protein (Cat# TP314598). The protein was produced from HEK293T cells transfected with NTHL1 cDNA clone (Cat# [RC214598]) using MegaTran 2.0 (Cat# [TT210002]).

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US