

Product datasheet for TP314651L

Neuroigin 4 (NLGN4Y) (NM_014893) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human neuroigin 4, Y-linked (NLGN4Y), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA	>RC214651 protein sequence
Clone or AA Sequence:	Red=Cloning site Green=Tags(s)

MLRPQGLLWLPPLLFTSVCVMLNSNVLLWITALAIKFTLIDSQAQYPVWNTNYGKIQGLRTPLPSEILGPV
EQYLGVPYASPTGERRFQPPESSWTGIRNATQFSAVCPQHLLDERFLLHDMPLPIWFTTSLDTLMTYVQ
DQNEDECLYLNIIYVPMEDDIHEQNSKKPVMVYIHGGSYMEGTGNMIDGSILASYGNVIVITINYRLGILGF
LSTGDQAAKGNYGLLDQIQALRWIENVGAFGGDPKRVTIFGSGAGASCVSLLTSLSHYSEGLFQKAIQS
GTALSSWAVNYQPAKYTRILADKVGCMMLDTTDMVECLKNKNYKELIQQTITPATYHIAFGPVIDGDVIP
DDPQILMEQGEFLNYDIMLGVNQGEGLKFVDGIVDNEDGVTPNDFDFSNSFVDNLYGYPEGKDTLRETI
KFMYTDWADKENPETRRKTLVALFDHQWVAPAVATADLHAQYGSPTYFYAFYHHCQSEMKPSWADSAHG
DEVYVFGIPMIGPTLFCNFSKNDVMSAVVMTYWTFNAKTGDPNQVVPQDTKFIHTKPNRFEEVAWS
KYNPKDQLYLHIGLKPRVRDHYRATKVAFWLELVPHLHNLNEIFQYVSTTTKVPPDMSFPYGTRRSPA
KIWPTTKRPAITPANNPKHSDPHKTGPEDTTVLIETKRDYSTELSVTIavgasllflnilafaalykk
DKRRHETHRHSPQRNTTNDITHIQNEEIMSLQMKQLEHDHECESLQAHDTLRLTCLPPDYTLRLRRSPDD
IPFMTPTNTITMIPNTLMGMQPLHTFKTFSGGQNSTNLPHGHSTTRV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	91.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



[View online »](#)

Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_055708
Locus ID:	22829
UniProt ID:	Q8NFZ3 , A0A024R9E7
RefSeq Size:	5409
Cytogenetics:	Yq11.221
RefSeq ORF:	2448
Synonyms:	HNL4Y
Summary:	This gene encodes a type I membrane protein that belongs to the family of neuroligins, which are cell adhesion molecules present at the postsynaptic side of the synapse, and may be essential for the formation of functional synapses. Alternatively spliced transcript variants have been found for this gene.[provided by RefSeq, Mar 2011]
Protein Families:	Druggable Genome, Transmembrane

Product images:



Coomassie blue staining of purified NLGN4Y protein (Cat# [TP314651]). The protein was produced from HEK293T cells transfected with NLGN4Y cDNA clone (Cat# [RC214651]) using MegaTran 2.0 (Cat# [TT210002]).