

Product datasheet for **TP314738M**

ZNF205 (NM_003456) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens zinc finger protein 205 (ZNF205), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC214738 protein sequence Red =Cloning site Green =Tags(s)
	<p>MSADGGGIQDTQDKETPPEVPDRGHPHQEMP SKLGEAVPSGDAQESLHIKMEPEEPHSE GASQEDGAQGA WGWAPLSHGSKEKALFLPGGALPSPRIPVLSREGRTRDRQMAAALLTAW SQMPVTFEDVALYLSREEWGR LDHTQQNFYRDVLQKKNGLSLGFPSRPFWAPQAHGKGEASGSSRQAGDEKEWRGACTGAVEVGQRVQTS SVAALGNVKPFRTRAGRVQWGV PQCAQEAACGRSSGPAKDSGQPAEPDRTPDAAPPDPSPTPEQEYRVPE KPNEEEEKGAPESGEEGLAPDSEVGRKSYRCEQCGKGFWSHSHLVTHRRHTHTGEKPYACTDCGKRFRSSH LIQHQQIHTGEKPYTCPACRKSFSHHSTLIQHQRHTHTGEKPYVCDRCAKRFTRRSDLVTHQGTHTGAKPH KCPICAKCFTQSSALVTHQRTHTGVKPYPCPEGKCFQSRSNLI AHNRTHTGEKPYHCLDCGKSFSSH LTAHQRTHRGVPRPYACPLCGKSFRRSNLHRHEKIHTTGPKALAMLMLGAAAAGALATPPPAPT</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	60.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_003447
Locus ID:	7755
UniProt ID:	O95201
RefSeq Size:	2077
Cytogenetics:	16p13.3
RefSeq ORF:	1662
Synonyms:	RhitH; Zfp13; ZNF210
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Protein Families:	Transcription Factors

Product images:



Coomassie blue staining of purified ZNF205 protein (Cat# [TP314738]). The protein was produced from HEK293T cells transfected with ZNF205 cDNA clone (Cat# [RC214738]) using MegaTran 2.0 (Cat# [TT210002]).