

## Product datasheet for **TP314894L**

### UBL7 (NM\_201265) Human Recombinant Protein

#### Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ubiquitin-like 7 (bone marrow stromal cell-derived) (UBL7), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC214894 representing NM_201265 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	<p>MSLSDWHLAVKLADQPLTPKSILRLPETELGEYSLGGYSISFLKQLIAGKQLQESVPDPELIDLIYCGRKL KDDQTLDFYGIQPGSTVHVLKRSWPEPDQKPEPVDKVAAMREFRVLHTALHSSSSYREAVFKMLSNKESL DQIIVATPGLSSDPIALGVLQDKDLFSVFADPNMLD TLVPAHPALVNAIVLVLH SVAGSAPMPGTDSSSR SMPSSSYRDMPGGFLFEGLSDEDDDFHPNTRSTPSSSTPSSRPASLGYSGAAGPRPITQSELATALALAS TPESSSHTPTPGTQGHSSGTSPMSSGVQSGTPITNDLFSQALQHALQASGQPSLQSQWQPQLQQLRDMGI QDDELSLRALQATGGDIQAALELIFAGGAP</p> <p><b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b></p>
Tag:	C-Myc/DDK
Predicted MW:	40.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_957717</a></u>



[View online »](#)

Locus ID:	84993
UniProt ID:	<a href="#">Q96S82</a>
RefSeq Size:	1409
Cytogenetics:	15q24.1
RefSeq ORF:	1140
Synonyms:	BMSC-UbP; TCBA1
Protein Families:	Druggable Genome

### Product images:



Coomassie blue staining of purified UBL7 protein (Cat# [TP314894]). The protein was produced from HEK293T cells transfected with UBL7 cDNA clone (Cat# [RC214894]) using MegaTran 2.0 (Cat# [TT210002]).