

## Product datasheet for TP314919

### Striatin 4 (STRN4) (NM\_001039877) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens striatin, calmodulin binding protein 4 (STRN4), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC214919 representing NM_001039877 Red=Cloning site Green=Tags(s)

MMEERAAAAVAAAASSCRPLGSGAGPGPTGAAPVSAPAPGPGPAGKGGGGGGSPGPTAGPEPLSLPGILH  
FIQHEWARFEAEKARWEAERAELQAQVAFLQGERKGQENLKTDLVRRIKMLEYALKQERAKYHKLKFGTD  
LNQGEKKADVSEQVSNPVSVTLENSPLVWKEGRQLLRQYLEEVGYTDTILDMRSKRVRSLGRSLELN  
GAVEPSEGAPRAPPGPAGLSGGESLLVKQIEEQIKRNAAGKDGERLGGSVLGQIPFLQNCEDSDSDDED  
ELDSVQHKKQRVKLPKALVPEMEDEDEEDDSEDAINEFDLFGSGEDGEGAPDRRCTVDGSPHELESRR  
VKLQGILADLRDVGDLPPKVTGPPPGTPQPRPHEGSFGFSSDVFIMDTIGGGEVSLGLDADLTVTNDNDL  
SCDLSDSKDAFKKTWNPKFTRLRSHYDGIRSLAFHHSQSALLTASEDGTLLKWNLQKAVTAKKNAALDVEP  
IHAFRAHRGPVLAVAMGSNSEYCYSGGADACIHSWKIPDLSMDPYDGYDPSVLSHVLEGHGDAVWGLAFS  
PTSQRLASCSADGTVRIWDPSSSSPACLCTFPTASEHGVPVTSVAFTSTEPAHIVASFRSGDVTLYDMEVG  
SALLTLESRGSSGPTQINQVSHPNQPLTITAHDDRIRFLDNRTGKPVHSMVAHLDAVTCLAVDPNGAF  
LMSGSHDCSLRLWSLDNKTCVQEITAHRRKKHEEAIHAVACHPSKALIASAGADALAKVFV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	81.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



[View online »](#)

Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_001034966</a>
Locus ID:	29888
UniProt ID:	<a href="#">Q9NRL3</a> , <a href="#">Q59GV6</a>
RefSeq Size:	3242
Cytogenetics:	19q13.32
RefSeq ORF:	2280
Synonyms:	PPP2R6C; ZIN; zinedin
Summary:	Binds calmodulin in a calcium dependent manner. May function as scaffolding or signaling protein.[UniProtKB/Swiss-Prot Function]

### Product images:



Coomassie blue staining of purified STRN4 protein (Cat# TP314919). The protein was produced from HEK293T cells transfected with STRN4 cDNA clone (Cat# [RC214919]) using MegaTran 2.0 (Cat# [TT210002]).