

Product datasheet for **TP315018L**

AGL (NM_000643) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins
Description: Recombinant protein of human amylo-1, 6-glucosidase, 4-alpha-glucanotransferase (AGL), transcript variant 3, 1 mg
Species: Human
Expression Host: HEK293T
Expression cDNA >RC215018 representing NM_000643
Clone or AA Sequence: Red=Cloning site Green=Tags(s)


MGHSKQIRILLNEMEKLEKTLFRLEQGYELQFRLGPTLQGKAVTVYTNYPFPGETFNREKFRSLDWENP
 TEREDDSKYCKLNLQSGSFQYYFLQGNEKSGGGYIVVDPILRVGADNHVLPDLCVTLQTFLLAKCLGPF
 DEWESRLRVAKESGYNMIHFTPLQTLGLSRSCYSLANQLELNPDFSRRPNRKYTWNDVGQLVEKLLKEWNV
 ICITDVVYNHTAANSKWIQEHPECAYNLVNSPHLKPWVLDLALWRFSCDVAEGKYKEKGIPALIENDHH
 MNSIRKIIWEDIFPKLKLWEFFQVDVKNKAVEQFRLLTQENRRVTKSDPNQHLTIIQDPEYRRFGCTVDM
 NIALTTFIPHDKGPAAIECCNWFHKRMEELNSEKHRLINYHQEQAVNCLLGNVYERLAGHGPKLGPVT
 RKHPLVTRYFTFPFEEIDFSMEESMIHLPNKACFLMAHNGWVMGDDPLRNFAEPGSEVYLRRELICWGDS
 VKLRYGNKPEDCPYLWAHMKKYTEITATYFQGVRLDNCHSTPLHVAEYMLDAARNLQPNLYVVAELFTGS
 EDLDNVFVTRLGISSLIREAMSAYNSHEEGRVYRYGGEPVGSFVQPCRPLMPAIAHALFMDITHDNEC
 PIVHRSAYDALPSTTIVSMACCASGSTRGYDELVPHQISVSEERFYTKWNPEALPSNTGEVNFQSGIIA
 ARCAISKLHQELGAKGFIQVYVDQVDEDIVAVTRHSPSIHQSVAVSRTAFRNPKTSFYSKEVPQMCIPG
 KIEEVLEARTIERNTKPYRKDENSINGTPDITVEIREHIQLNESKIVKQAGVATKGPNEYIQEIEFENL
 SPGSVIIFRVSLDPAHQVAVGILRNHLTQFSPHFKSGSLAVDNADPILKIPFASLASRLTLAELNQILYR
 CESEEKEDGGGCYDIPNWSALKYAGLQGLMSVLAIRPKNDLGHFPFCNNLRSGDWMIDYVSNRLISRSGT
 IAEVGKWLQAMFFYLKQIPRYLIPCYFDAILIGAYTLLDTAWKQMSSFVQNGSTFVKHLSLGSVQLCGV
 GKFPSLPILSPALMDVPYRLNEITKEKEQCCVSLAAGLPHFSSGIFRCWGRDFTIALRGILLITGRYVEA
 RNIILAFAGTLRHGLIPNLLGEGIYARYNCRDAVWWWLQCIQDYCKMVPNGLDILKCPVSRMYPTDDSDAP
 LPAGTLDQPLFEVIQEAMQKHMQGIQFRERNAGPQIDRNMKDEGFNITAGVDEETGFVYGGNRFNCGTWM
 DKMGESDRARNRGIPATPRDGSVEIVGLSKSAVRWLELSKKNIFPYHEVTVKRHGKAIKVSYDEWNRK
 IQDNFEKLFHVSEDPDLNEKHPNLVHKRGIYKDSYGASSPWCDYQLRPNFTIAMVVAPELFTTEKAWKA
 LEIAEKLLGPLGMKTLDPDDMVYCGIYDNALDNDNYNLAKGFNYHQPEWLWPIGYFLRAKLYFSRLMG
 PETTAKTIVLKNVLSRHYVHLERSPWKGLPELTNENAQYCPFCSETQAWSIATILETLTYDL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online >](#)

Tag:	C-Myc/DDK
Predicted MW:	174.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_000634
Locus ID:	178
UniProt ID:	P35573 , A0A0S2A4E4
RefSeq Size:	7165
Cytogenetics:	1p21.2
RefSeq ORF:	4596
Synonyms:	GDE
Summary:	This gene encodes the glycogen debrancher enzyme which is involved in glycogen degradation. This enzyme has two independent catalytic activities which occur at different sites on the protein: a 4-alpha-glucotransferase activity and a amylo-1,6-glucosidase activity. Mutations in this gene are associated with glycogen storage disease although a wide range of enzymatic and clinical variability occurs which may be due to tissue-specific alternative splicing. Alternatively spliced transcripts encoding different isoforms have been described. [provided by RefSeq, Jul 2008]
Protein Families:	Druggable Genome
Protein Pathways:	Metabolic pathways, Starch and sucrose metabolism

Product images:

Coomassie blue staining of purified AGL protein (Cat# [TP315018]). The protein was produced from HEK293T cells transfected with AGL cDNA clone (Cat# [RC215018]) using MegaTran 2.0 (Cat# [TT210002]).