

Product datasheet for **TP315035L**

DNAJC12 (NM_201262) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human Dnaj (Hsp40) homolog, subfamily C, member 12 (DNAJC12), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC215035 representing NM_201262 Red =Cloning site Green =Tags(s)

MDAILNYRSEDTEDYYTLLGCDELSSVEQILAEFKVRALECHPDKHPENPKAVETFQKLQKAKEILTNEE
SRARYDHWRRSQMSMPFQQWEALNDSVKTVGFSLGAT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	12.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_957714
Locus ID:	56521
UniProt ID:	Q9UKB3
RefSeq Size:	786



[View online »](#)


Cytogenetics: 10q21.3

RefSeq ORF: 321

Synonyms: HPANBH4; JDP1

Summary: This gene encodes a member of a subclass of the HSP40/Dnaj protein family. Members of this family of proteins are associated with complex assembly, protein folding, and export. Two transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]

Product images:



Coomassie blue staining of purified DNAJC12 protein (Cat# [TP315035]). The protein was produced from HEK293T cells transfected with DNAJC12 cDNA clone (Cat# [RC215035]) using MegaTran 2.0 (Cat# [TT210002]).