

Product datasheet for TP315181

PAX3 (NM_181457) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens paired box 3 (PAX3), transcript variant PAX3, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC215181 representing NM_181457 Red =Cloning site Green =Tags(s) MTTLAGAVPRMMRPGPGQNYPRSGFPLEVSTPLGQGRVNQLGGVFINGRPLPNHIRHKIVEMAHHGIRP C VISRQLRVSHGCVSKILCRYQETGSIRPGAIGGSKPKQVTTDPVEKKIEEYKRENPGMFSWEIRDKLLKD AVCDRNTVPSVSSISRILRSKFGKGEEEEADLERKEAESEKKAKHSIDGILSERASAPQSDEGSDIDSE PDLPLKRKQRRSRTTFTAEQLEELERAFERTHYPDITYREELAQRAKLTEARVQVWFNRRARWRKQAGA NQLMAFNHLIPGGFPPTAMPTLPTYQLSETSYQPTSIPQAVSDPSSTVHRPQLPPSTVHQSTIPSNPDS SSAYCLPSTRHGFSSYTDSFVPPSGPSNPMNPTIGNGLSPQVMGLLTNHGGVPHQPQTDYALSPLTGGLE PTTTVSASCQRLDHMKSLDSLPTSQSYCPPTYSTTGYSMDPVTGYQYGYGQSKPWTF TR TRPLE QKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	52.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.


[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_852122</u>
Locus ID:	5077
UniProt ID:	<u>P23760</u>
RefSeq Size:	3193
Cytogenetics:	2q36.1
RefSeq ORF:	1437
Synonyms:	CDHS; HUP2; WS1; WS3
Summary:	This gene is a member of the paired box (PAX) family of transcription factors. Members of the PAX family typically contain a paired box domain and a paired-type homeodomain. These genes play critical roles during fetal development. Mutations in paired box gene 3 are associated with Waardenburg syndrome, craniofacial-deafness-hand syndrome, and alveolar rhabdomyosarcoma. The translocation t(2;13)(q35;q14), which represents a fusion between PAX3 and the forkhead gene, is a frequent finding in alveolar rhabdomyosarcoma. Alternative splicing results in transcripts encoding isoforms with different C-termini. [provided by RefSeq, Jul 2008]
Protein Families:	Adult stem cells, Druggable Genome, Embryonic stem cells, ES Cell Differentiation/IPS, Transcription Factors

Product images:



Coomassie blue staining of purified PAX3 protein (Cat# TP315181). The protein was produced from HEK293T cells transfected with PAX3 cDNA clone (Cat# [RC215181]) using MegaTran 2.0 (Cat# [TT210002]).