

Product datasheet for **TP315182**

FAM151B (NM_205548) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human family with sequence similarity 151, member B (FAM151B), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC215182 protein sequence Red =Cloning site Green =Tags(s)

MAASAGGPGSWSENILEYFLRNSQITAEDGAEITWYHAANHKAQTNEALKSTAHMIEADVLLPSDGEHS
QPIMAHPPETNSDNTLQEWLTEVMKSNKGIKLDKSLAVVEPSMMLLENVKRHLKRPVWINADILPGPNG
NSKVIDAKPFLDVISFFPDVTFSLGWTTGWHPEKVNESWTMVKEMEYICNELSQPVTFPVRAALVRQ
SCSQLLWLLKSNRYSLTIWTGKNDNYSVEDLLYIRDHFDDKKQVFYDILEPQNHEFKQAIGIKVNL


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	31.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_991111</u>
Locus ID:	167555



[View online »](#)

UniProt ID: [Q6UXP7](#)
RefSeq Size: 1707
Cytogenetics: 5q14.1
RefSeq ORF: 828
Synonyms: UNQ9217

Product images:

Coomassie blue staining of purified FAM151B protein (Cat# TP315182). The protein was produced from HEK293T cells transfected with FAM151B cDNA clone (Cat# [RC215182]) using MegaTran 2.0 (Cat# [TT210002]).