

## Product datasheet for **TP315203M**

### ETS1 (NM\_005238) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human v-ets erythroblastosis virus E26 oncogene homolog 1 (avian) (ETS1), transcript variant 2, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC215203 representing NM_005238 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MKAAVDLKPTLTIKTEKVDLELFPSPDMECADVPLLTPSSKEMMSQALKATFSGFTKEQQRLGIPKDPR  
QWTETHVRDWMWAVNEFSLKGVDFQKFCMNGAALCALGKDCFLELAPDFVGDILWEHLEILQKEDVKPY  
QVNGVNPAYPESRYTSDYFISYGIEHAQCVPPSEFSEPSFITESYQTLHPISSEELLSLKYENDYPSVIL  
RDPLQDTLQNDYFAIKQEVTPDNMCMGRTSRGLGGQDSFESIESYDSCDRLTQSWSSQSSFNSLQRV  
PSYDSFDSYDPAALPNHKPKGTFKDYVRDRADLNKDKPVIPAAALAGYTGSGPIQLWQFLLELLTDKSC  
QSFISWTGDGWFKLSDPDEVARRWGRKRNKPKMNYEKLSRGLRYYYDKNIIHKTAGKRYVYRFVCDLQS  
LLGYTPEELHAMLVDKPDADA

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**



Tag:	C-Myc/DDK
Predicted MW:	50.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol



[View online »](#)

<b>Bioactivity:</b>	ETS1 Activity Verified in a DNA-binding Assay: <b>ETS1 (TP315203)</b> activity was measured in a colorimetric DNA-binding assay. Purified <b>ETS1</b> protein containing a C-terminal MYC/DDK tag was incubated with biotinylated double-stranded oligonucleotide containing the ETS1 consensus DNA-binding sequence. Following incubation, the reaction was transferred to a streptavidin-coated microplate to allow capture of the protein-DNA complex. After washing, the captured protein was detected with an anti-DDK peroxidase conjugate and colorimetric signal detection with TMB. Specificity of the protein-DNA interaction was confirmed by carrying out the binding in the presence of an unlabeled competitor oligonucleotide and by comparison to binding to an oligonucleotide containing a mutation in the consensus binding sequence.
<b>Preparation:</b>	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
<b>Note:</b>	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
<b>Storage:</b>	Store at -80°C.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_005229</a>
<b>Locus ID:</b>	2113
<b>UniProt ID:</b>	<a href="#">P14921</a> , <a href="#">B4DW78</a>
<b>RefSeq Size:</b>	5228
<b>Cytogenetics:</b>	11q24.3
<b>RefSeq ORF:</b>	1323
<b>Synonyms:</b>	c-ets-1; ETS-1; EWSR2; p54
<b>Summary:</b>	This gene encodes a member of the ETS family of transcription factors, which are defined by the presence of a conserved ETS DNA-binding domain that recognizes the core consensus DNA sequence GGAA/T in target genes. These proteins function either as transcriptional activators or repressors of numerous genes, and are involved in stem cell development, cell senescence and death, and tumorigenesis. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.[provided by RefSeq, Jul 2011]
<b>Protein Families:</b>	Druggable Genome, Transcription Factors
<b>Protein Pathways:</b>	Dorso-ventral axis formation, Pathways in cancer, Renal cell carcinoma

Product images:



Coomassie blue staining of purified ETS1 protein (Cat# [TP315203]). The protein was produced from HEK293T cells transfected with ETS1 cDNA clone (Cat# [RC215203]) using MegaTran 2.0 (Cat# [TT210002]).