

## Product datasheet for TP315233L

### HDAC7 (NM\_015401) Human Recombinant Protein

#### Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human histone deacetylase 7 (HDAC7), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC215233 representing NM_015401 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MHSPGADGTQVSPGAHYCSPTGAGCPRPCADTPGPPQPMPDLRVGQRPPVEPPPPEPTLLALQRPQRLHHH LFLAGLQQRSVPEMRLSMDTPMPELQVGPQEQELRQLLHKDKSKRSVASSVVKQKLAEVILKKQQAAL ERTVHPNSPGIPYRTLEPLETEGATRSMLSSFLPPVPSLPSDPPEHFPLRKTVSEPNLKLRYKPKKSLER RKNPLLRKESAPPSLRRRPAETLGDSSPSSSTPASGCSSPNDSEHGPNPILGSEALLGQRLRLQETSVA PFALPTVSLLPAILTGLPAPARADSDRRTHPTLGRGPILGSPHTPLFLPHGLEPEAGGTLPSRLQPILL LDPSGSHAPLLTVPGLGPLPFHFAQSLMTTERLSGSGLHWPLSRTRSEPLPPSATAPPPGPMQPRLEQL KTHVQVIKRSKAKPSEKPLRQIPSAEDLETGGGPGQVDDGLEHRELGHGQPEARGPAPLQQHPQVLLW EQQLLAGRLPRGSTGDTVLLPLAQGGHRPLSRAQSSPAAPASLSAEPASQARVLSSETPARTLPFTTG LIYDSVMLKHQCSCGDNSRHPEHAGRIQSIWSRLQERGLRSQCECLRGRKASLEELQSVHSERHVLLYGT NPLSRLKLDNGKLAGLLAQRMFVMLPCGGVGVDTDTIWNELHSSNAARWAAGSVTDLAFKVASRELKNGF AVRPPGHHADHSTAMGFCCFNVAIACRQLQQSKASKILIVDWDVHHGNGTQQTFYQDPSVLYISLHR HDDGNFFPGSGAVDEVGAGSGEGFNVNVAWAGGLDPPMGDPEYLAFRIVMPIAREFSPDLVLVSAGFD AAEGHPAPLGGYHVSACFCGYMTQQLMNLGAGVLALEGGHDLTAICDASEACVAALLGNRVDPLSEEG WKQKPNLNAIRSLEAVIRVHSKYWGCMQRLASCPDSWVPRVPGADKEEVEAVTALASLSVGILAEDRPSE QLVEEEPMNL
	<b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	106.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol



[View online »](#)

<b>Preparation:</b>	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
<b>Note:</b>	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
<b>Storage:</b>	Store at -80°C.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<u>NP_056216</u>
<b>Locus ID:</b>	51564
<b>UniProt ID:</b>	<u>Q8WUI4</u>
<b>RefSeq Size:</b>	4301
<b>Cytogenetics:</b>	12q13.11
<b>RefSeq ORF:</b>	2973
<b>Synonyms:</b>	HD7; HD7A; HDAC7A
<b>Summary:</b>	Histones play a critical role in transcriptional regulation, cell cycle progression, and developmental events. Histone acetylation/deacetylation alters chromosome structure and affects transcription factor access to DNA. The protein encoded by this gene has sequence homology to members of the histone deacetylase family. This gene is orthologous to mouse HDAC7 gene whose protein promotes repression mediated via the transcriptional corepressor SMRT. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]
<b>Protein Families:</b>	Druggable Genome, Transcription Factors

### Product images:



Coomassie blue staining of purified HDAC7 protein (Cat# [TP315233]). The protein was produced from HEK293T cells transfected with HDAC7 cDNA clone (Cat# [RC215233]) using MegaTran 2.0 (Cat# [TT210002]).