

Product datasheet for TP315335

Protein Kinase D2 (PRKD2) (NM_016457) Human Recombinant Protein

Product data:

| | |
|---------------------------------------|---|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant protein of human protein kinase D2 (PRKD2), transcript variant 1, 20 µg |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | >RC215335 representing NM_016457 Red=Cloning site Green=Tags(s) |

MATAPSYAGLPGSPGPGSPPPPGGLELQSPPLLPQIPAPGSGVSFHIQIGLTREFVLLPAASELAHVK
QLACSIVDQKFPECGFYGLYDKILLFKHDPTSANLLQLVRSSGDIQEGDLVEVWLSASATFEDFQIRPHA
LTVHSYRAPAFCDHCGEMLFGLVRQGLKCDGCLNYHKRCAFSIPNNCSGARKRRLSSTSLASGHSVRLG
TSESLPCTAEELSRSTTELLPRRPPSSSSSSASSYTRPIELDKMLLSKVKVPHTFLIHSYTRPTVCQA
CKKLLKGLFRQGLQCKDCKFNCHKRCATRVPNDCLEALINGDVPMEEATDFSEADKSALMDESEDSGVI
PGSHSENALHASEEEEEGEGGKAQSSLGYIPLMRVVQSVRHTRKSTTLREGWVHYSNKDTLRKRHYWR
LDCKITLQFNNTTNRYYKEIPLSEILTVESAQNFLVPPGTNPHCFEIVTANATYFVGEMPPGGTGGPS
GQGAEAARGWETAIRQALMPVILQDAPSAPGHAPHRQASLSISVSNSQIQENVDIATVYQIFPDEVLGSG
QFGVVYGGKHKRKTGRDVAVKVIDKLRFPKQESQLRNEVAILQSLRHPGIVNLECMFETPEKVFVMEKL
HGDMLEMILSSEKGRLEPERLTKFLITQILVALRHLHFKNIVHCDLKPENVLLASADPPQVKLCDFGFAR
IIGKSFRRSVVGTAPYLAPEVLLNQGYNRSLDMWSVGVIMYVSLSGTFPFNEDEDINDQIQNAAFMYPA
SPWSHISAGAILINLLQVKMRKRYSDKSLSHPWLQEYQTWLDLRELEGKMGERYITHESDDARWEQF
AAEHPLPGSGLPTDRDLGGACPPQDHMQGLAERISVL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

| | |
|----------------|--|
| Tag: | C-Myc/DDK |
| Predicted MW: | 96.5 kDa |
| Concentration: | >0.05 µg/µL as determined by microplate BCA method |
| Purity: | > 80% as determined by SDS-PAGE and Coomassie blue staining |
| Buffer: | 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol |
| Preparation: | Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps. |



[View online »](#)

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_057541](#)

Locus ID: 25865

UniProt ID: [Q9BZL6](#)

RefSeq Size: 2883

Cytogenetics: 19q13.32

RefSeq ORF: 2634

Synonyms: HSPC187; nPKC-D2; PKD2

Summary: The protein encoded by this gene belongs to the protein kinase D (PKD) family of serine/threonine protein kinases. This kinase can be activated by phorbol esters as well as by gastrin via the cholecystokinin B receptor (CCKBR) in gastric cancer cells. It can bind to diacylglycerol (DAG) in the trans-Golgi network (TGN) and may regulate basolateral membrane protein exit from TGN. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome, Protein Kinase

Product images:



Coomassie blue staining of purified PRKD2 protein (Cat# TP315335). The protein was produced from HEK293T cells transfected with PRKD2 cDNA clone (Cat# [RC215335]) using MegaTran 2.0 (Cat# [TT210002]).