

Product datasheet for **TP315353L**

FAM72D (NM_207418) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human gastric cancer up-regulated-2 (GCUD2), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC215353 protein sequence Red =Cloning site Green =Tags(s)
	MSTNICSFKDRCVSILCCKFCKQVLSSRGMKAVLLADTEIDLFSTDIPPTNAVDFTGRCYFTKICKCKLK DIACLKCGNIVGYHVIVPCSSCLLSCNNGHFWMFHSQAVYDINRLDSTGVNVLLWGNLPEIEESTDEDVL NISAEECIR
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	16.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_997301</u>
Locus ID:	728833
UniProt ID:	<u>Q6L9T8</u>
RefSeq Size:	2393



[View online »](#)

Cytogenetics: 1q21.1

RefSeq ORF: 447

Synonyms: GCUD2

Product images:



Coomassie blue staining of purified FAM72D protein (Cat# [TP315353]). The protein was produced from HEK293T cells transfected with FAM72D cDNA clone (Cat# [RC215353]) using MegaTran 2.0 (Cat# [TT210002]).