

Product datasheet for TP315387

NVL (NM_002533) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human nuclear VCP-like (NVL), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC215387 representing NM_002533 Red=Cloning site Green=Tags(s)

MKPRPAGFVDNKLKQRVIQYLTSNCKGKYVDIGVLASDLQRVYSIDYGRRKRNAFRIQVEKVFSIISSEK
ELKNLLEDEHLAKRARQGEEDNEYTESYSDDDSSMEDYDPDPQSANHMNSSLLSLYRKGNPDSVSNTPE
MEQRETTSSSTPRISSKTGSIPLKTPAKDSEGGWFIDKTPSVKKDSFFLDLSCEKSNPKKPIEQDSKDS
SLESMDKRKGLKNKGSKRKEDLQEVDGEIEAVLQKKAKARGLEFQISNVKFEDVGGNDMTLKEVCKM
LIHMRHPEVYHHLGVPPRGVLLHGPPGCGKTLAHAIAGELDLPILKVAAPEIVSGVSGESEQKLRELF
EQAVSNAPCIIFIDEIDAITPKREVASKDMERRIVAQLLTCMDDLNNVAATARVLVIGATNRPDSDLPAL
RRAGRFDREICLGIPDEASRERILQTLCRKLRLPQAFDFCHLAHLTPGFVGADLMALCREAAMCAVNRVL
MKLQEQQKNPEMEDLPKGVQEERLGTEPTSETQDELQRLGLLRDQDPLSEEQMQGLCIELNDFIVAL
SSVQPSAKREGFVTPNVTWADIGALEDIREELTMAILAPVRNPDQFKALGLVTPAGVLLAGPPGCGKTL
LAKAVANESGLNFISVKPELLNMYVGESERAVRQVFQRAKNSAPCVIFFDEVDALCPRRSRETGASVR
VVNQLLTEMDGLEARQQVFIMAATNRPDIIDPAILRPGRLDKTLFVGLPPPADRLAILKTITKNGTKPPL
DADVNL EAIAGDLRCDCTGADLSALVREASICALRQEMARQKSGNEKGELKVSHKHFEFAFKKVRSSIS
KKDQIMYERLQESLSR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	94.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_002524](#)

Locus ID: 4931

UniProt ID: [O15381](#)

RefSeq Size: 2921


Cytogenetics: 1q42.11

RefSeq ORF: 2568

Synonyms: NVL2

Summary: This gene encodes a member of the AAA (ATPases associated with diverse cellular activities) superfamily. Multiple transcript variants encoding different isoforms have been found for this gene. Two encoded proteins, described as major and minor isoforms, have been localized to distinct regions of the nucleus. The largest encoded protein (major isoform) has been localized to the nucleolus and shown to participate in ribosome biosynthesis (PMID: 15469983, 16782053), while the minor isoform has been localized to the nucleoplasm. [provided by RefSeq, Aug 2011]

Product images:



Coomassie blue staining of purified NVL protein (Cat# TP315387). The protein was produced from HEK293T cells transfected with NVL cDNA clone (Cat# [RC215387]) using MegaTran 2.0 (Cat# [TT210002]).