

Product datasheet for **TP315402**

FUBP3 (NM_003934) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens far upstream element (FUSE) binding protein 3 (FUBP3), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC215402 representing NM_003934 Red =Cloning site Green =Tags(s)
	<p>MAELVQGSAPVGMKAEGFVDALHRVRQIAAKIDSIPHLNNSPLVDPSVYGVYQKRPLDDGVDGNQLGALVHQRVTITEEFKVPDKMVGFIIGRGGEQISRIQAESGCKIQIASESSGIPERPCVLGTGTPESIEQAKRL LGQIVDRCRNGPGFHNDIDSNSTIQEILIPASKVGLVIGRGGETIKQLQERTGVKMVMIQDGPLPTGADK PLRITGDAFKVQQAREMVLEIIREKDQADFRGVRGDFNSRMGGGSIEVSVPRFAVGIVIGRNGEMIKKIQ NDAGVRIQFKPDDGISPERAAQVMGPPDRCQAAHIISELILTAQERDGFGLAAARGRGRGRGDWSVGA PGGVQEITYTVPADKCGLVIGKGGENIKSINQQSGAHVELQRNPPNSDPNLRRTIRGVPQQIEVARQL IDEKVGGTNLGAPGAFGQSPFSQPPAPPHQNTFPPRSSGCFPNMAAKVNGNPHSTPVSQPAPFLTQGWGS TYQAWQQPTQQVPSQQSQPQSSQPNYSKAWEDYYKKQSHAASAAPQASSPPDYTMAWAEYYRQQVAFYGGQ TLGQAQAHSQEQ</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	61.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_003925
Locus ID:	8939
UniProt ID:	Q96I24 , A0A024R8A7
RefSeq Size:	3175
Cytogenetics:	9q34.11-q34.12
RefSeq ORF:	1716
Synonyms:	FBP3
Summary:	May interact with single-stranded DNA from the far-upstream element (FUSE). May activate gene expression.[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified FUBP3 protein (Cat# TP315402). The protein was produced from HEK293T cells transfected with FUBP3 cDNA clone (Cat# [RC215402]) using MegaTran 2.0 (Cat# [TT210002]).