

Product datasheet for TP315525

HYAL1 (NM_007312) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human hyaluronoglucosaminidase 1 (HYAL1), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC215525 representing NM_007312 Red =Cloning site Green =Tags(s)
	<p>MAAHLIPICALFLTLDDMAQGFRGPLLPNRPFTTVWNANTQWCLERHGVDDVDSVFDVVANPGQTFRGPD MTIFYSSQLGTYPYTPTGEPVFGGLPQNASLIAHLARTFQDILAAIPAPDFSGLAVIDWEAWRPRWAFN WDTKDIYRQRSRALVQAQHPDWPAPQVEAVAQDQFQGAARAWMAGTLQLGRALRPRGLWGFYGFPCYNY DFLSPNYTGQCPSGIRAQNDQLGWLWGQSRALYPSIYMPAVLEGTGKSQMYVQHRVAEAFRVAAGA DPNLPVLPYVQIFYDTTNHFLPLDELEHSLGESAAQGAAGVVLWVSWENTRTKESCQAIKEYMDTTLGPFILN VTSGALLCSQALCSGHGRCVRRTSHPKALLLNPAFSIQLTPGGGPLSLRGALSLEDQAQMAVEFKCRC YPGWQAPWCERKSMW</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	48.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_009296



[View online »](#)

Locus ID: 3373

UniProt ID: [Q12794](#)

RefSeq Size: 2518

Cytogenetics: 3p21.31

RefSeq ORF: 1305


Synonyms: HYAL-1; LUCA1; MGC45987; NAT6

Summary: This gene encodes a lysosomal hyaluronidase. Hyaluronidases intracellularly degrade hyaluronan, one of the major glycosaminoglycans of the extracellular matrix. Hyaluronan is thought to be involved in cell proliferation, migration and differentiation. This enzyme is active at an acidic pH and is the major hyaluronidase in plasma. Mutations in this gene are associated with mucopolysaccharidosis type IX, or hyaluronidase deficiency. The gene is one of several related genes in a region of chromosome 3p21.3 associated with tumor suppression. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

Protein Families: Secreted Protein

Protein Pathways: Glycosaminoglycan degradation, Lysosome, Metabolic pathways

Product images:



Coomassie blue staining of purified HYAL1 protein (Cat# TP315525). The protein was produced from HEK293T cells transfected with HYAL1 cDNA clone (Cat# [RC215525]) using MegaTran 2.0 (Cat# [TT210002]).