

Product datasheet for TP315549L

GRIN3A (NM_133445) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human glutamate receptor, ionotropic, N-methyl-D-aspartate 3A (GRIN3A), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC215549 representing NM_133445 Red=Cloning site Green=Tags(s)

MRRLSLWWLLSRVCLLLPPPCALVLAGVPSSSSHPQPCQILKRIGHAVRVGAVHLQPWTTAPRAASRAPD
DSRAGAQRDEPEPGTRRSPAPSPGARWLGSTLHGRGPPGSRKPGEGARAEALWPRDALLFAVDNLNRVEG
LLPYNLSLEVVMAIEAGLGDPLLPFSSPSPWSSDPFSFLQSVCHTVVQGVSAALLAFPQSQGEMMELD
LVSLVLHIPVISIVRHEFPRESQNPPLHLQLSLENSLSSDADVTVSILTMNNWYNFSLLLCQEDWNITDFL
LLTQNNKFHLGSIINITANLPSTQDLLSFLQIQLESIKNSTPTVVMFGCDMESIRRIFEITTQFGVMPP
ELRWVLGDSQNVEELRTEGLPLGLIAHGKTTQSVFEHYVQDAMELVARAVATATMIQPELALIPSTMNCM
EVETTNTLSGQYLSRFLANTTFRGLSGSIRVKGSTIVSSENNFFIWNLQHDPMGKPMWTRLGSWQGGKIV
MDYGIWPEQAQRHKTHFQHPSKLHLRVVTLIEHPFVFTREVDDEGLCPAGQLCLDPMTNDSSTLDSLFS
LHSSNDTVPIKFKCCYGYCIDLLEKIAEDMNFDFDLYIVGDGKYGAWKNGHWTLVGDLLRGTAHMAVT
SFSINTARSQVIDFTSPFFSTSLGILVRTRDTAAPIGAFMWPLHWTMWLGIFVALHITAVFLTYEWKSP
FGLTPKGRNRSKVFSSALNICYALLFGRTVAIKPPKCWTGRFLMNLWAIFCMFCLSTYTANLAAVMVG
EKIYEELSGIHDPKLHHPQSQGFRFGTVRESSAEDYVRQSFPEMHEYMRRYNVPATPDGVEYLKNDPEKLD
AFIMDKALLDYEVSIDADCKLLTVGKPFIAIEGYGIGLPPNSPLTANISELISQYKSHGFMDMLHDKWYRV
VPCGKRFAVTETLQMGIKHFSGLFVLLCIGFGLSILTTIGEIVYRLLLPRIKNSKLQYWLHTSQRLH
RAINTSFIEEKQQHFKTKRVEKRSNVGPRQLTVWNTSNLSDNRRKYIFSDEEQNQLGIQIHQDIPLPP
RRRELPAALRTTNGKADSLNVSRSVMQELSELEKQIQVIRQELQLAVSRKTELEEYQRTSRTCES

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	125.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol



[View online »](#)

Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_597702
Locus ID:	116443
UniProt ID:	Q8TCU5
RefSeq Size:	7770
Cytogenetics:	9q31.1
RefSeq ORF:	3345
Synonyms:	GluN3A; NMDAR-L; NMDAR3A; NR3A
Summary:	This gene encodes a subunit of the N-methyl-D-aspartate (NMDA) receptors, which belong to the superfamily of glutamate-regulated ion channels, and function in physiological and pathological processes in the central nervous system. This subunit shows greater than 90% identity to the corresponding subunit in rat. Studies in the knockout mouse deficient in this subunit suggest that this gene may be involved in the development of synaptic elements by modulating NMDA receptor activity. [provided by RefSeq, Jul 2008]
Protein Families:	Druggable Genome, Ion Channels: Glutamate Receptors, Transmembrane
Protein Pathways:	Neuroactive ligand-receptor interaction

Product images:



Coomassie blue staining of purified GRIN3A protein (Cat# [TP315549]). The protein was produced from HEK293T cells transfected with GRIN3A cDNA clone (Cat# [RC215549]) using MegaTran 2.0 (Cat# [TT210002]).