

Product datasheet for **TP315631L**

Ecat1 (KHDC3L) (NM_001017361) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 6 open reading frame 221 (C6orf221), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC215631 representing NM_001017361 Red =Cloning site Green =Tags(s)
	MDAPRRFPTLVQLMQPKAMPVEVLGHLPKRFSWFHSEFLKNPKVVRLEVWLVEKIFGRGGERIPHVQGM S QILIHVNRLDPNGEAEILVFGRPYQEDTIKMIMNLADYHRQLQAKGSGKALAQDVATQKAETQRSSIEV REAGTQRSVEVREAGTQRSVEVQEVGTQGSPVEVQEAGTQQSLQAANKSGTQRSPEAASKAVTQRFREDA RDPVTRL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	24.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001017361</u>
Locus ID:	154288



[View online »](#)

UniProt ID: [Q587J8](#)

RefSeq Size: 1063


Cytogenetics: 6q13

RefSeq ORF: 651

Synonyms: C6orf221; ECAT1; HYDM2

Summary: The protein encoded by this gene belongs to the KHDC1 family, members of which contain an atypical KH domain that may not bind RNA like canonical KH domains. This gene is specifically expressed in the oocytes, and recent studies suggest that it may function as a regulator of genomic imprinting in the oocyte. Mutations in this gene are associated with recurrent biparental complete hydatidiform mole. [provided by RefSeq, Dec 2011]

Product images:



Coomassie blue staining of purified KHDC3L protein (Cat# [TP315631]). The protein was produced from HEK293T cells transfected with KHDC3L cDNA clone (Cat# [RC215631]) using MegaTran 2.0 (Cat# [TT210002]).