

Product datasheet for **TP315636**

MACROD2 (NM_001033087) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human MACRO domain containing 2 (MACROD2), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC215636 protein sequence Red =Cloning site Green =Tags(s)
	<p>MNEFFSVDDNNEEEEDVEMKEDSDENGPEEKQSVEEMEEQSQDADGVNTVTVPGPASEEAVEDCKDEDFA KDENITKGGVTDHSVRDQDHPDGQENDSTKNEIKIETESQSSYMETEELSSNQEDAVIVEQPEVIPLTE DQEEKEGEGAPGEDTPRMPGKSEGSSDLENTPGPDVEMNSQVDKVNTPTESQQEDQLIAGAQDEAKEQRN GTK</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	23.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001028259
Locus ID:	140733



[View online »](#)

UniProt ID: [A1Z1Q3](#)

RefSeq Size: 4425


Cytogenetics: 20p12.1

RefSeq ORF: 639

Synonyms: C2orf133; C20orf133

Summary: The protein encoded by this gene is a deacetylase involved in removing ADP-ribose from mono-ADP-ribosylated proteins. The encoded protein has been shown to translocate from the nucleus to the cytoplasm upon DNA damage. [provided by RefSeq, May 2017]

Product images:



Coomassie blue staining of purified MACROD2 protein (Cat# TP315636). The protein was produced from HEK293T cells transfected with MACROD2 cDNA clone (Cat# [RC215636]) using MegaTran 2.0 (Cat# [TT210002]).