

Product datasheet for TP315649

SUNC1 (SUN3) (NM_001030019) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human Sad1 and UNC84 domain containing 1 (SUNC1), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC215649 representing NM_001030019 Red =Cloning site Green =Tags(s)
	 MSGKTKARRAAMFFRRCSEDASGSASGNALLSEDENPDANGVTRSWKIILSTMLTLTFLLVGLLNHQWLK ETDVPQKSRQLYAIIEYGSRLYKYQARLRMPKEQLELLKKESQNLNFRQILFLIEQIDVLKALLRDM KDGMDNNHNWNTHGDPVEDPDHTEEVSNLVNYVLKLRDQVEMADYALKSAGASIIIEAGTSESYKNNKA KLYWHGIGFLNHEMPPDIILQPDVYPGKCWAFPGSQGHTLIKLATKIIPAVTMEHISEKVSPSGNISSA PKFESVYGITKKCEGEEIFLGQFIYNKTGTTVQTFELQHAVSEYLLCVKLNIFSNWGHGPKYTCLYRFRVH GTPGKHI TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	40.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001025190



[View online »](#)

Locus ID: 256979

UniProt ID: [Q8TAQ9](#)

RefSeq Size: 1452


Cytogenetics: 7p12.3

RefSeq ORF: 1071

Synonyms: SUNC1

Summary: As a probable component of the LINC (LInker of Nucleoskeleton and Cytoskeleton) complex, involved in the connection between the nuclear lamina and the cytoskeleton. The nucleocytoplasmic interactions established by the LINC complex play an important role in the transmission of mechanical forces across the nuclear envelope and in nuclear movement and positioning. May be involved in nuclear remodeling during sperm head formation in spermatogenesis. A probable SUN3:SYNE1 LINC complex may tether spermatid nuclei to posterior cytoskeletal structures such as the manchette.[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified SUN3 protein (Cat# TP315649). The protein was produced from HEK293T cells transfected with SUN3 cDNA clone (Cat# [RC215649]) using MegaTran 2.0 (Cat# [TT210002]).