

Product datasheet for TP315772L

BCAS4 (NM_017843) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human breast carcinoma amplified sequence 4 (BCAS4), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC215772 representing NM_017843 Red =Cloning site Green =Tags(s)
	<p>MQRTGGGAPRPRGNHGLPGSLRQPDVALLMLLVADADQPEPMRSGARELALFLTPEPGAEAKEVEETIEG MLLRLEEFCSLADLIRSQTILEENIPVLKAKLTEMRGYAKVDRLEAFVKMVGHHVAFLEADVLAQER DHGAFPQALRRWLGSAGLPSFRNVECSGTIPARCNLRLPGSSDSPASASQVAGITEVTCTGARDVRAAHT V</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	22.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_060313</u>
Locus ID:	55653



[View online »](#)

UniProt ID: [Q8TDM0](#)
RefSeq Size: 1368
Cytogenetics: 20q13.13
RefSeq ORF: 633
Synonyms: CNOL

Product images:

Coomassie blue staining of purified BCAS4 protein (Cat# [TP315772]). The protein was produced from HEK293T cells transfected with BCAS4 cDNA clone (Cat# [RC215772]) using MegaTran 2.0 (Cat# [TT210002]).