

Product datasheet for **TP315831M**

HIST1H2AH (NM_080596) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human histone cluster 1, H2ah (HIST1H2AH), 100 µg

Species: Human

Expression Host: HEK293T

**Expression cDNA Clone
or AA Sequence:** >RC215831 protein sequence
Red=Cloning site **Green**=Tags(s)

MSGRGKQGGKARAKAKTRSSRAGLQFPVGRVHLLRKGNYAERVGAGAPVYLAHVLEYLTAEILELAGNA
ARDNKKTRIIPRHLQLAIRNDEELNKLKLVKVTIAQGGVLPNIQAVLLPKKTESHHKAK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 13.7 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_542163](#)

Locus ID: 85235

UniProt ID: [Q96KK5](#), [A3KPC7](#)

RefSeq Size: 481

Cytogenetics: 6p22.1



[View online »](#)


RefSeq ORF: 384

Synonyms: dj86C11.1; H2A/S; H2AFALii; H2AH; HIST1H2AH

Summary: Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. Two molecules of each of the four core histones (H2A, H2B, H3, and H4) form an octamer, around which approximately 146 bp of DNA is wrapped in repeating units, called nucleosomes. The linker histone, H1, interacts with linker DNA between nucleosomes and functions in the compaction of chromatin into higher order structures. This gene is intronless and encodes a replication-dependent histone that is a member of the histone H2A family. Transcripts from this gene lack polyA tails but instead contain a palindromic termination element. This gene is found in the histone microcluster on chromosome 6p21.33. [provided by RefSeq, Aug 2015]

Protein Pathways: Systemic lupus erythematosus

Product images:



Coomassie blue staining of purified HIST1H2AH protein (Cat# [TP315831]). The protein was produced from HEK293T cells transfected with HIST1H2AH cDNA clone (Cat# [RC215831]) using MegaTran 2.0 (Cat# [TT210002]).