

Product datasheet for **TP315869M**

KLHL30 (NM_198582) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human kelch-like 30 (Drosophila) (KLHL30), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC215869 representing NM_198582 Red =Cloning site Green =Tags(s)

MVRNVDDLDFHLP SHAQDMLDGLQRLRSQPKLADVTLLVGGRELPC HRGLLALSSPYFHAMFAGDFAE
SF
SARVELRDVEPAVVGQLVDFVYTGRLTITQGNVEALTRTAARLHFPSVQKVCGRYLQQQLDAANCLGICE
FGEQQGLLGVA AKAWAFLRENFEAVAREDEFLQLPRERLVTCLAGDLLQVQPEQSRLEALMRWVRHDPQ
A
RAAHLPELLSLVHLDVAPRPCVQQLLASEPLIQESEACRAALSQGH DGAPLALQQKLEEVLVVVG GQALE
EEEAGEEPTPGLGNFAFYNSKAKRWMALPDFPDYHKWGFSLAALNNNIYVTGGSRGKTDTWSTTQAWC
F
PLKEASWKPVAPMLKPRTNHAS AALNGEIYVIGGTTLDVVEVESYDPYTD SWTPVSPALKYVS NFSAAGC
RGRLYLVGSSACKYNALALQCYNPVTDAWSVIASPFLPKYLSSPRCAALHGELYLIGDNTKKVYVYDPGA
NLWQKVQS QHSLHENGALVPLGDALYVTGGRWQGM EGDYHVEMEAYDTRD TWTRHGALPRLWLYH
GAST
VFLDVSKWTQPSGPTQEH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	63.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_940984](#)

Locus ID: 377007


UniProt ID: [Q0D2K2](#)

RefSeq Size: 3557

Cytogenetics: 2q37.3

RefSeq ORF: 1734

Product images:



Coomassie blue staining of purified KLHL30 protein (Cat# [TP315869]). The protein was produced from HEK293T cells transfected with KLHL30 cDNA clone (Cat# [RC215869]) using MegaTran 2.0 (Cat# [TT210002]).