

Product datasheet for **TP315910L**

MAGIX (NM_024859) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human MAGI family member, X-linked (MAGIX), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC215910 representing NM_024859 Red =Cloning site Green =Tags(s)

MEPRTGGAANPKGSRGSRGPSPLAGPSARQLLARLDARPLAARAAMDVAALVRRAGATLRLRRKEAVSVL
DSADIEVTD SRLPHATIVDHRPQHRWLETNAPPQLIQGKARSAPKPSQASGHFSVELVRGYAGFGLTLG
GGRDVAGDTPLAVRGLLKDGPAQRRCGRLEVGDVLVHINGESTQGLTHAQAVERIRAGGPQLHLVIRRRPLE
THPGKPRGVGEPKRGVWPDRSPDPGGPEVTGSRSSSTSLVQHPPSRITLKKTRGSPEPSPEAAADGP
TVSPERRAEDPNDQIPGSPGPWLVPSEERLSRALGVRGAAQLAQEMAAGRRRH


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	35 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_079135
Locus ID:	79917



[View online »](#)

UniProt ID: [Q9H6Y5](#)
RefSeq Size: 2150
Cytogenetics: Xp11.23
RefSeq ORF: 1002
Synonyms: JM10; PDZX

Product images:

Coomassie blue staining of purified MAGIX protein (Cat# [TP315910]). The protein was produced from HEK293T cells transfected with MAGIX cDNA clone (Cat# [RC215910]) using MegaTran 2.0 (Cat# [TT210002]).