

Product datasheet for TP315972

RBM47 (NM_019027) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human RNA binding motif protein 47 (RBM47), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC215972 representing NM_019027 Red =Cloning site Green =Tags(s)
	<p>MTAEDSTAAMSSDSAAGSSAKVPEGVAGAPNEAALLALMERTGYSMVQENGQRKYGGPPPGWEGPHPQRG CEVFGKIPRDVYEDELVPVFEAVGRIYELRLMMDFDGKNRGYAFVMYCHKHEAKRAVRELNYYEIRPGR LLGVCCSVDNCRLLFIGGIPKMKKREEILEEIAKVTEGVLDVIVYASAADKMKNRGFAFVEYESHRAAAMA RRKLMPGRIQLWGHQIAVDWAEPEIDVDEDVDMETVKILYVRNLMIETTEDTIKKSFGQFNPGCVERVKKI RDYAFVHFTSREDAVHAMNNLNGTELEGSCLEVTLAKPVDKEQYSRYQKAARGGGAAEAAQQPSVYVYSCD PYTLAYYGYPYNALIGPNRDYFVKVAIPAIGAQYSMFPAAPAPKMIEDGKIHTVEHMISPIAVQPD PASA AAAAAAAAAAAAAVIPTVSTPPPFQGRPITPVYTVAPNVQRIPTAGIYGASYVPPAAPATATIATLQKNA AAAAAVYGGYAGYIPQAFPAAAIQVPIPDVYQTY</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	56.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP_061900](#)

Locus ID: 54502

UniProt ID: [A0AV96](#)


RefSeq Size: 4020

Cytogenetics: 4p14

RefSeq ORF: 1572

Synonyms: NET18

Product images:



Coomassie blue staining of purified RBM47 protein (Cat# TP315972). The protein was produced from HEK293T cells transfected with RBM47 cDNA clone (Cat# [RC215972]) using MegaTran 2.0 (Cat# [TT210002]).