

Product datasheet for TP316129

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

HYAL1 (NM_033159) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human hyaluronoglucosaminidase 1 (HYAL1), transcript variant 7, 20 μg

Species: Human
Expression Host: HEK293T

Expression cDNA Clone or AA Sequence: >RC216129 protein sequence Red=Cloning site Green=Tags(s)

 $MAAHLLPICALFLTLLDMAQGFRGPLLPNRPFTTVWNANTQWCLERHGVDVDVSVFDVVANPGQTFRGPD\\ MTIFYSSQLGTYPYYTPTGEPVFGGLPQNASLIAHLARTFQDILAAIPAPDFSGLAVIDWEAWRPRWAFN$

WDTKDIYRQRSRALVQAQHPDWPAPQVEAVAQDQFQGAARAWMAGTLQLGRALRPRGLWGFYGFPDCYNY DFLSPNYTGQCPSGIRAQNDQLGWLWGQSRALYPSIYMPAVLEGTGKSQMYVQHRVAEAFRVAVAAGDPN LPVLPYVQIFYDTTNHFLPLDELEHSLGESAAQGAAGVVLWVSWENTRTKESCQAIKEYMDTTLGPFILN VTSGALLCSQALCSGHGRCVRRTSHPKALLLLNPASFSIQLTPGGGPLSLRGALSLEDQAQMAVEFKCRC

YPGWQAPWCERKSMW

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK
Predicted MW: 48.2 kDa

Concentration: >0.05 μg/μL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional

chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some

loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling

conditions. Avoid repeated freeze-thaw cycles.

RefSeq: NP 149349





Locus ID: 3373

UniProt ID: Q12794, A0A024R2X3

RefSeq Size: 2103 Cytogenetics: 3p21.31 RefSeq ORF: 1305

Synonyms: HYAL-1; LUCA1; MPS9; NAT6

Summary: This gene encodes a lysosomal hyaluronidase. Hyaluronidases intracellularly degrade hyaluronan,

one of the major glycosaminoglycans of the extracellular matrix. Hyaluronan is thought to be involved in cell proliferation, migration and differentiation. This enzyme is active at an acidic pH


and is the major hyaluronidase in plasma. Mutations in this gene are associated with

mucopolysaccharidosis type IX, or hyaluronidase deficiency. The gene is one of several related genes in a region of chromosome 3p21.3 associated with tumor suppression. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

Protein Families: Secreted Protein

Protein Pathways: Glycosaminoglycan degradation, Lysosome, Metabolic pathways

Product images:



Coomassie blue staining of purified HYAL1 protein (Cat# TP316129). The protein was produced from HEK293T cells transfected with HYAL1 cDNA clone (Cat# [RC216129]) using

MegaTran 2.0 (Cat# [TT210002]).