

Product datasheet for **TP316151L**

Von Hippel Lindau (VHL) (NM_000551) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human von Hippel-Lindau tumor suppressor (VHL), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC216151 representing NM_000551 Red =Cloning site Green =Tags(s)
	MPRRAENWDEAEVGAEEAGVEEYGPPEEDGGEEESGAEEESGPPEESGPEELGAEEMEAGRPRPVLRVNSRE PSQVIFCNRSRVLVPWLNFDGEPQPYPTLPPGTGRRHSYRGHLWLF RDAGTHDGLLVNQTELFVPSL NVDGQPIFANITLPVYTLKERCLQVVRSLVKPENYRRLDIVRSLYEDLEDHPNVQKDLERLTQERIAHQR MGD
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	24 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Bioactivity:	Binding assay (PMID: 27780863)
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_000542



[View online »](#)


Locus ID: 7428
UniProt ID: [P40337](#)
RefSeq Size: 2968
Cytogenetics: 3p25.3
RefSeq ORF: 639
Synonyms: HRCA1; pVHL; RCA1; VHL1

Summary: Von Hippel-Lindau syndrome (VHL) is a dominantly inherited familial cancer syndrome predisposing to a variety of malignant and benign tumors. A germline mutation of this gene is the basis of familial inheritance of VHL syndrome. The protein encoded by this gene is a component of the protein complex that includes elongin B, elongin C, and cullin-2, and possesses ubiquitin ligase E3 activity. This protein is involved in the ubiquitination and degradation of hypoxia-inducible-factor (HIF), which is a transcription factor that plays a central role in the regulation of gene expression by oxygen. RNA polymerase II subunit POLR2G/RPB7 is also reported to be a target of this protein. Alternatively spliced transcript variants encoding distinct isoforms have been observed. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome, Transcription Factors

Protein Pathways: Pathways in cancer, Renal cell carcinoma, Ubiquitin mediated proteolysis

Product images:



Coomassie blue staining of purified VHL protein (Cat# [TP316151]). The protein was produced from HEK293T cells transfected with VHL cDNA clone (Cat# [RC216151]) using MegaTran 2.0 (Cat# [TT210002]).