

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# Product datasheet for TP316151M

### Von Hippel Lindau (VHL) (NM\_000551) Human Recombinant Protein

#### **Product data:**

Product Type:	Recombinant Proteins		
Description:	Recombinant protein of human von Hippel-Lindau tumor suppressor (VHL), transcript variant 1, 100 μg		
Species:	Human		
Expression Host:	HEK293T		
Expression cDNA Clone or AA Sequence:	>RC216151 representing NM_000551 <mark>Red</mark> =Cloning site Green=Tags(s)		
	MPRRAENWDEAEVGAEEAGVEEYGPEEDGGEESGAEESGPEESGPEELGAEEEMEAGRPRPVLRSVNSRE PSQVIFCNRSPRVVLPVWLNFDGEPQPYPTLPPGTGRRIHSYRGHLWLFRDAGTHDGLLVNQTELFVPSL NVDGQPIFANITLPVYTLKERCLQVVRSLVKPENYRRLDIVRSLYEDLEDHPNVQKDLERLTQERIAHQR MGD		
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV		
Tag:	C-Myc/DDK		
Predicted MW:	24 kDa		
Concentration:	>0.05 µg/µL as determined by microplate BCA method		
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining		
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol		
Bioactivity:	Binding assay (PMID: <u>27780863</u> )		
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.		
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.		
Storage:	Store at -80°C.		
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.		
RefSeq:	<u>NP 000542</u>		



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	Von Hippel Lindau (VHL) (NM_000551) Human Recombinant Protein – TP316151M	
Locus ID:	7428	
UniProt ID:	<u>P40337</u> , <u>A0A024R2F2</u>	
RefSeq Size:	2968	
Cytogenetics:	3p25.3	
RefSeq ORF:	639	
Synonyms:	HRCA1; pVHL; RCA1; VHL1	
Summary:	Von Hippel-Lindau syndrome (VHL) is a dominantly inherited familial cancer syndrome predisposing to a variety of malignant and benign tumors. A germline mutation of this gene is the basis of familial inheritance of VHL syndrome. The protein encoded by this gene is a component of the protein complex that includes elongin B, elongin C, and cullin-2, and possesses ubiquitin ligase E3 activity. This protein is involved in the ubiquitination and degradation of hypoxia-inducible-factor (HIF), which is a transcription factor that plays a central role in the regulation of gene expression by oxygen. RNA polymerase II subunit POLR2G/RPB7 is also reported to be a target of this protein. Alternatively spliced transcript variants encoding distinct isoforms have been observed. [provided by RefSeq, Jul 2008]	
Protein Families:	Druggable Genome, Transcription Factors	
Protein Pathways	S: Pathways in cancer, Renal cell carcinoma, Ubiquitin mediated proteolysis	

## **Product images:**

116 —	-
66 —	
45 —	H
35 —	-
25 —	
18 —	
14	

Coomassie blue staining of purified VHL protein (Cat# [TP316151]). The protein was produced from HEK293T cells transfected with VHL cDNA clone (Cat# [RC216151]) using MegaTran 2.0 (Cat# [TT210002]).

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US