

Product datasheet for **TP316161L**

PRX (NM_020956) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human periaxin (PRX), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC216161 representing NM_020956 Red =Cloning site Green =Tags(s) MEARSRSAEELRRAELVEIIVETEAQTGVSGINVAGGGKEGIFVRELREDS PAARSLSLQEGDQLLSARV FFENFKYEDALRLLQCAEPYKVSFCLKRTVPTGDLALRPGTVSGYEIKGPRAKVAKL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	15.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_066007
Locus ID:	57716
UniProt ID:	Q9BXM0
RefSeq Size:	5502
Cytogenetics:	19q13.2



[View online »](#)


RefSeq ORF: 383

Synonyms: CMT4F

Summary: This gene encodes a protein involved in peripheral nerve myelin upkeep. The encoded protein contains 2 PDZ domains which were named after PSD95 (post synaptic density protein), DlgA (Drosophila disc large tumor suppressor), and ZO1 (a mammalian tight junction protein). Two alternatively spliced transcript variants have been described for this gene which encode different protein isoforms and which are targeted differently in the Schwann cell. Mutations in this gene cause Charcot-Marie-Tooth neuropathy, type 4F and Dejerine-Sottas neuropathy. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified PRX protein (Cat# [TP316161]). The protein was produced from HEK293T cells transfected with PRX cDNA clone (Cat# [RC216161]) using MegaTran 2.0 (Cat# [TT210002]).