

Product datasheet for **TP316229L**

NHLRC2 (NM_198514) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human NHL repeat containing 2 (NHLRC2), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC216229 representing NM_198514 Red =Cloning site Green =Tags(s)

MAAPGGRGRSLSGLLPAQTSLEYALLDAVTQQEKDSLVIYQYLQKVDGWEQDLSVPEFPEGLEWLNTEEPI
SVYKDLGKIVLDFFTYCCINCIIHLLPDLHALEHTYSKDKGLLIIGVHSAKFPNEKVLNLIKSAVLRYN
ITHPMVNDADASLWQELEVSCWPTLVILGPRGNMLFSLIGEGHKDKLFLYTSIALKYYKDRGQIRDNKIG
IKLYKDSLPPSPLLPFGKVTVDQVTDRLVIADTGHHRILVWVKNGQIQYSIGGPNPGRKDGIFSESTFNS
PQGVAIMNNIIYVADTENHLIRKIDLEAEKVSTVAGIGIQGTDKEGGAKGEQQPISSPWDVVFVGTSGSEV
QRGDILWIAMAGTHQIWALLLDGKLPKKNELTKGTCLRFAGSGNEENRNNAYPHKAGFAQPSGLSLASE
DPWSCLVADSESSTVRTVSLKDGAVKHLVGGGERDPMNLFAGDVGVDGGINAKLQHPLGVTWDKKRNLly
VADSYNHKIKVDPKTKNCTTLAGTGDTNNVTSSSFTESTFNEPGGLCIGENGELLYVADTNNHQIKVMD
LETKMVSVLPPIRSENAVVDGPFLVEKQKTLPKLPKSAPSIRLSPVTACAGQTLQFKLRLDLPVSGSKLTE
GVSSCWFLTAEGNEWLLQGQIAAGDIENISSQPTISLQIPDDCLSLEAIVSVSVFLYCSADSSACMMKA
ILFSQPLQITDTQQGCIAPVELRYVF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	79.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_940916
Locus ID:	374354
UniProt ID:	Q8NBF2 , Q7Z658
RefSeq Size:	2550
Cytogenetics:	10q25.3
RefSeq ORF:	2178
Synonyms:	FINCA
Summary:	Required for normal embryonic development.[UniProtKB/Swiss-Prot Function]
Protein Families:	Druggable Genome

Product images:



Coomassie blue staining of purified NHLRC2 protein (Cat# [TP316229]). The protein was produced from HEK293T cells transfected with NHLRC2 cDNA clone (Cat# [RC216229]) using MegaTran 2.0 (Cat# [TT210002]).