

Product datasheet for **TP316422L**

PITPNM2 (NM_020845) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human phosphatidylinositol transfer protein, membrane-associated 2 (PITPNM2), 1 mg
Species:	Human
Expression Host:	HEK293T



[View online »](#)

Expression cDNA Clone or AA Sequence: >RC216422 representing NM_020845
Red=Cloning site Green=Tags(s)


MIIKEYRIPLPMTVEEYRIAQLYMIQKKSERNETYGEGSGVEILENRPYTDGPGGSGQYTHKVYHVGMHIP
 SWFRSILPKAALRWVEESWNAYPYTRTRFTCPFVEKFSIDIETFYKTDAGENPDVFNLSPEKNQLTIDF
 IDIVKDPVPHNEYKTEEDPKLFQSTKTQRGPLSENWIEEYKKQVFPIPCAYKLCKVEFRYWGMQSKIERF
 IHDTGLRRVMVRAHRQAWCWQDEWYGLSMENIRELEKEAQLMLSRKMAQFNEDGEEATELVKHEAVSD
 QT
 SGEPPEPSSSNGEPLVGRGLKKQWSTSSKSSSRSSKRGASPSRHSISEWRMQSIARDSDESSDDEFFDAHE
 DLSLTEEMFPKDITKWSSNDLMDKIESPEPEDTQDGLYRQGAPEFRVASSVEQLNIIEDVVSQPLAAPP
 KIHVLLLVHGGTILDTGAGDPSSKKGDANTIANVFDTVMRVHYPALGRLAIRLVPCPPVCSDAFALVS
 NLSPYSHDEGLSSSQDHIPLAALPLLATSSPQYQEAATVIQRANLAYGDFIKSQEGMTFNGQVCLIGD
 CVGGILAFDALCYSNQPVSSESQSSSRGSSVSMQDNDLLSPGILMNAAHCCGGGGGGGGGGSSGGGG
 SS
 GGSSLESSRHLRSNVDIPRSNGTEDPKRQLPRKRSDSSTYELDTIQQHQAFLSSLHASVLRTEPCSRHS
 SSSTMLDGTGALGRFDFEITDLFLFGCPLGLVLALRKTVPALDVFQLRPACQQVYNLFHPADPSASRL
 PLLERRFHALPPFSVPRYQRYPLGDGCSTLLADVLQTHNAAFQEHGAPSSPGTAPASRGFRRASEISIAS
 QVSGMAESYTASSIAQKAPDALSHTPSVRRLLSLLALPAPSPTTPGPHPPARKASPGLERAPGLPELDIGE
 VAAKWWGQKRIDYALYCPDALTAFTVALPHLFHASYWESTDVVSFLLRQVMRHDNSSILELDGKEVSFV
 TPSKPREKWQRKRTHVKLRNVTANHRINDALANEDGPQVLTGRFMYGPLDMVTLTGEKVDVHIMTQPP
 SG
 EWLYLDTLVTNNSGRVSYTIPESHRLGVGVYPIKMVVRGDHTFADSYITVLPKGTEFVVSIDGSFAASV
 SIMGSDPKVRAGAVDVVRHWQDLGYLIIVTGRPDMQKQRRVVAWLAQHNFPHGVVSFCDGLVHDLPLRH
 KA
 NFLKLLISELHLRVHAAYGSTKDVAVYSAISLSPMQIYVGRPTKKLQQQCQFITDGYAAHLAQLKYSHR
 ARPARNATRMALRKGSFGLPGQGDFLRSRNHLLRTISAQPSGSPSHRHERTQSQADGEQGRQSRMSVAA
 G
 CWGRAMTGRLEPGAAAGPK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK
Predicted MW: 148.8 kDa
Concentration: >0.05 µg/µL as determined by microplate BCA method
Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage: Store at -80°C.

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_065896</u>
Locus ID:	57605
UniProt ID:	<u>Q9BZ72</u>
RefSeq Size:	6754
Cytogenetics:	12q24.31
RefSeq ORF:	4047
Synonyms:	NIR-3; NIR3; RDGB2; RDGBA2
Summary:	PITPNM2 belongs to a family of membrane-associated phosphatidylinositol transfer domain-containing proteins that share homology with the Drosophila retinal degeneration B (rdgB) protein (Ocaka et al., 2005 [PubMed 15627748]).[supplied by OMIM, Mar 2008]
Protein Families:	Druggable Genome

Product images:



Coomassie blue staining of purified PITPNM2 protein (Cat# [TP316422]). The protein was produced from HEK293T cells transfected with PITPNM2 cDNA clone (Cat# [RC216422]) using MegaTran 2.0 (Cat# [TT210002]).