

Product datasheet for TP316631

C14orf184 (NM_001080113) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human hypothetical protein LOC650662 (PP8961), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC216631 representing NM_001080113 Red =Cloning site Green =Tags(s)
	MSVPKLVNHGWETPCNPDPALLLTGAVPGEGETTLERALPRVVGEQELGSWRERANSTPERRLSPALKH VSARPTSIPETPGERGARRRDRAGRARKGEGAAASTLSSLRGAAPAEARGAEAPESVSGGADAVPLEG APAVRSKS
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	15.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001073582
Locus ID:	650662
RefSeq Size:	2596
Cytogenetics:	14q32.12




[View online »](#)

RefSeq ORF: 444

Synonyms: FLJ33842; PP8961

Summary: DISCONTINUED: This record has been withdrawn by NCBI staff. This record was defined by AF318368.1 and AK091161.1, which appear to be cloning artifacts and are not sufficient evidence to define a distinct locus.

Product images:



Coomassie blue staining of purified C14orf184 protein (Cat# TP316631). The protein was produced from HEK293T cells transfected with C14orf184 cDNA clone (Cat# [RC216631]) using MegaTran 2.0 (Cat# [TT210002]).