

## Product datasheet for **TP316704M**

### **MNK2 (MKNK2) (NM\_199054) Human Recombinant Protein**

#### Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human MAP kinase interacting serine/threonine kinase 2 (MKNK2), transcript variant 2, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC216704 representing NM_199054 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	<p>MVQKKPAELQGFHRSFKGQNPFEAFSLDQPDHGDSDFLQCSARPDMPASQPIDIPDAKKRGKKKKRGR</p> <p>ATDSFSGRFEDVYQLQEDVLGEGA HARVQTCINLITSQEYAVKIEKQPGHIRSRVFREVMELYQCQGHR</p> <p>NVLELIEFFEEEDRFYLVFEKMRGGSILSHIHKRRHFNELEASVVQDVASALDFLHNKGIAHRDLKPEN</p> <p>ILCEHPNQVSPVKICDFDLGSGIKLNGDCSPISTPELLTPCGSAEYMAPEVVEAFSEEASIYDKRCDLWS</p> <p>LGVILYILLSGYPPFVGRCGSDCGWDRGEACPACQNMLFESI QEGKYEPDKDWAHISCAAKDLISKLLV</p> <p>RDAKQRLSAAQVLQHPWVQGCAPENTLPTPMVLQRNSCAKDLTSFAAEAIAMNRQLAQHDEDLAEFEA</p> <p>AG</p> <p>QGQPVLVRATSRCLQLSPPSQSKLAQRRQRASLSSAPVVLVDHA</p> <p><b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b></p>
Tag:	C-Myc/DDK
Predicted MW:	51.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<u>NP_951009</u>
<b>Locus ID:</b>	2872
<b>UniProt ID:</b>	<u>Q9HBH9</u>
<b>RefSeq Size:</b>	3795
<b>Cytogenetics:</b>	19p13.3
<b>RefSeq ORF:</b>	1395
<b>Synonyms:</b>	GPRK7; MNK2
<b>Summary:</b>	This gene encodes a member of the calcium/calmodulin-dependent protein kinases (CAMK) Ser/Thr protein kinase family, which belongs to the protein kinase superfamily. This protein contains conserved DLG (asp-leu-gly) and ENIL (glu-asn-ile-leu) motifs, and an N-terminal polybasic region which binds importin A and the translation factor scaffold protein eukaryotic initiation factor 4G (eIF4G). This protein is one of the downstream kinases activated by mitogen-activated protein (MAP) kinases. It phosphorylates the eukaryotic initiation factor 4E (eIF4E), thus playing important roles in the initiation of mRNA translation, oncogenic transformation and malignant cell proliferation. In addition to eIF4E, this protein also interacts with von Hippel-Lindau tumor suppressor (VHL), ring-box 1 (Rbx1) and Cullin2 (Cul2), which are all components of the CBC(VHL) ubiquitin ligase E3 complex. Multiple alternatively spliced transcript variants have been found, but the full-length nature and biological activity of only two variants are determined. These two variants encode distinct isoforms which differ in activity and regulation, and in subcellular localization. [provided by RefSeq, Aug 2011]
<b>Protein Families:</b>	Druggable Genome, Protein Kinase
<b>Protein Pathways:</b>	Insulin signaling pathway, MAPK signaling pathway

### Product images:



Coomassie blue staining of purified MNK2 protein (Cat# [TP316704]). The protein was produced from HEK293T cells transfected with MNK2 cDNA clone (Cat# [RC216704]) using MegaTran 2.0 (Cat# [TT210002]).