

Product datasheet for TP316724

HIF1 beta (ARNT) (NM_001668) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human aryl hydrocarbon receptor nuclear translocator (ARNT), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC216724 representing NM_001668 Red =Cloning site Green =Tags(s)
	<p>MAATTANPEMTSDVPSLGPAIASGNSGPGIQGGGAIVQRAIKRRPGLDFDDDGEGNSKFLRCDDDDQMSND KERFARSDDEQSSADKERLARENHSEIERRRRNKMTAYITELSDMVPTCSALARKPKDLTILRMAVSHMK SLRGTGNTSTDGSYKPSFLTDQELKHLILEAADGFLFVSCETGRVWVSDSVTPVLNQPQSEWFGSTLY DQVHPDDVDKLRQLSTSENALTGRILDLTGTVKKEGQQSSMRMCMGSRRSFICRMRCGSSSVDPVSVN RLSFVRNRCRNLGSKVDGEPHFVWHCTGYIKAWPPAGVSLPDDDDPEAGQGSKFCLVAIGRLQVTSSPN CTDMSNVCQPTEFISRHNIEGIFTVDHRCVATVGYQPQELLGKNIVEFCHPEDQQLLRDSFQQVWKLKG QVLSVMFRFRSKNQEWLWMRTSSFTFQNPYSDEIEYIICTNTNVKNSSQEPRTLNTIQRPQLGPTANL PLEMGSGQLAPRQQQQTELDMVPGRDGLASYNHSQVQPVTTTGPEHSPLEKSDGLFAQDRDPRFSEI YHNINADQSKGISSSTVPATQQLFSGNTFPPTPRPAENFRNSGLAPPVTIVQPSASAGQMLAQISRHSN PTQGATPTWPTTRSGFSAQQVATQATAKTRTSQFGVGSFQTPSSFSSMSLPGAPTASPGAAAAYPSLTNR GSNFAPETGQTAGQFQTRTAEGVGVWPQWQGGQPHHRSSSEQHVVQPPAQQPGQPEVFQEMLSMLGDQS NSYNNEEFDLTMFPPFSE</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	86.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001669
Locus ID:	405
UniProt ID:	P27540
RefSeq Size:	4846
Cytogenetics:	1q21.3
RefSeq ORF:	2367
Synonyms:	bHLHe2; HIF-1-beta; HIF-1beta; HIF1-beta; HIF1B; HIF1BETA; TANGO
Summary:	This gene encodes a protein containing a basic helix-loop-helix domain and two characteristic PAS domains along with a PAC domain. The encoded protein binds to ligand-bound aryl hydrocarbon receptor and aids in the movement of this complex to the nucleus, where it promotes the expression of genes involved in xenobiotic metabolism. This protein is also a co-factor for transcriptional regulation by hypoxia-inducible factor 1. Chromosomal translocation of this locus with the ETV6 (ets variant 6) gene on chromosome 12 have been described in leukemias. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2013]
Protein Families:	Druggable Genome, Transcription Factors
Protein Pathways:	Pathways in cancer, Renal cell carcinoma

Product images:



Coomassie blue staining of purified ARNT protein (Cat# TP316724). The protein was produced from HEK293T cells transfected with ARNT cDNA clone (Cat# [RC216724]) using MegaTran 2.0 (Cat# [TT210002]).