

Product datasheet for TP316773L

SH3RF3 (NM_001099289) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human SH3 domain containing ring finger 3 (SH3RF3), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC216773 representing NM_001099289 Red=Cloning site Green=Tags(s)

MLLGASWLCASKAAAAAAQSEGDEDRPGERRRRRAAATAAGAGEDMDESSLDDLLECSVCLERLDTTAKV
LPCQHTFCRRCLLESIVCSRHELRCPECRILVGCQVDEL PANILLVRLLDGIRQRPRAGTSPGGSPPARPI
PGQSAAPTLAGGGGGAAGSTPGSPVFLSAAAGSTAGSLRELATSRTAPAAKNPCLLPYGKALYSYEGKEP
GDLKFNKGDIIVLRRKVDEQWYHGELHGTQGF LPASIQICIQLPHAPPQ GKALYDFEMKDKDQDKDCLT
FTKDEILTVLRRVDENWAE GMLGDKIGIFPLLYVELNDSAKQLIEMDKPCPAAAASSCNASLPDSGAVAS
VAPSPTLSSSGAVSAFQRRVDGKKNTKKRHSFTALSVTHRSSQAASHRHSMEISAPVLISSSDPRAAARI
GDLAHLSCAAPTQDVSSSAGSTPTAVPRAASVS GEQGT PPKVQLPLNVYLALYAYKPQKSDELELHKGEM
YRVLEKCQDGFK GASLRTGVSGVFPGNVYTPVSRVPAGGAGPPRNNVVGGSPLAKGITTTMHPGSGSLS
SLATATRPALPITTPQAHAQHPTASPTGSLRHS AQPTASQARSTISTA AHSAAQAQDRPTATVSPLRT
QNSPSRLPATSLRPHSVVSPQ HSHQPPVQMCP RPAIPLTSAASAITPPNVSAANLNGEAGGGPIGVLSTS
SPTNTGCKLDEKKSEKKEKSGLLKLLAGASTKKKSRSPSVSPHDPQVAVDALLQGAVGPEVSSLSIH
GRAGSCPIESEMQGAMGMEPLHRKAGSLDLNFTSPSRQAPLSMAAIRPEPKLLPRERYRVVVSYPQSEA
EIELKEGDIVFVHKKREDGWYKGT LQRNGRTGLFPGSFVESF


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	92.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001092759
Locus ID:	344558
UniProt ID:	Q8TEJ3
RefSeq Size:	5608
Cytogenetics:	2q13
RefSeq ORF:	2646
Synonyms:	POSH2; SH3MD4
Summary:	Has E3 ubiquitin-protein ligase activity.[UniProtKB/Swiss-Prot Function]

Product images:

Coomassie blue staining of purified SH3RF3 protein (Cat# [TP316773]). The protein was produced from HEK293T cells transfected with SH3RF3 cDNA clone (Cat# [RC216773]) using MegaTran 2.0 (Cat# [TT210002]).