

## Product datasheet for TP316914L

### TSH beta (TSHB) (NM\_000549) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens thyroid stimulating hormone, beta (TSHB), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC216914 representing NM_000549 Red=Cloning site Green=Tags(s)
	MTALFLMSMLFGLACGQAMSFCIPTEYTMHIERRECAVCLTINTTICAGYCMTRDINGKLFLPKYALSQD VCTYRDFIYRTVEIPGCPLHVAPYFSYPVALSCKCGKCNLDYSDCIHEAIKTNCTKPQKSYLVGFSV
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	13.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_000540</a>
Locus ID:	7252
UniProt ID:	<a href="#">P01222</a>
RefSeq Size:	578



[View online »](#)

Cytogenetics: 1p13.2

RefSeq ORF: 414


Synonyms: TSH-B; TSH-BETA

**Summary:** The four human glycoprotein hormones chorionic gonadotropin (CG), luteinizing hormone (LH), follicle stimulating hormone (FSH), and thyroid stimulating hormone (TSH) are dimers consisting of alpha and beta subunits that are associated noncovalently. The alpha subunits of these hormones are identical, however, their beta chains are unique and confer biological specificity. Thyroid stimulating hormone functions in the control of thyroid structure and metabolism. The protein encoded by this gene is the beta subunit of thyroid stimulating hormone. Mutations in this gene are associated with congenital central and secondary hypothyroidism and Hashimoto's thyroiditis. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, May 2013]

**Protein Families:** Druggable Genome, Secreted Protein

**Protein Pathways:** Autoimmune thyroid disease, Neuroactive ligand-receptor interaction

### Product images:



Coomassie blue staining of purified TSHB protein (Cat# [TP316914]). The protein was produced from HEK293T cells transfected with TSHB cDNA clone (Cat# [RC216914]) using MegaTran 2.0 (Cat# [TT210002]).