

Product datasheet for TP317053L

DGKD (NM_152879) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins
Description: Recombinant protein of human diacylglycerol kinase, delta 130kDa (DGKD), transcript variant 2, 1 mg
Species: Human
Expression Host: HEK293T
Expression cDNA Clone or AA Sequence: >RC217053 representing NM_152879
Red=Cloning site Green=Tags(s)

MAAAAGAPPPGPPQPPPPPPPEESSDSEPEAEPGSPQKLIRKVSTSGQIRQKTIIEGMLTKQNNFQRS
 KRRYFKLRGRTLYYAKTAKSIIFDEVLDLDASVAESSTKNVNSFTVITPCRKLILCADNRKEMEDWIAA
 LKTVQNREHFPTQYSMDHFSGMHNHWYACSHARPTYCNVCREALSGVTS HGLSCEVCKFKAHKRAVRAT
 NNCKWTTLASIGKDIIEDADGIAMPHQWLEGNLPVSAKCTVCDKTCGSLVRLQLDWRCLWCKAMVHTSCKE
 SLLTKCPLGLCKVSVIPPTALNSIDSDGFWKASCPPSCTSPLLVFNKSGDNQGVKFLRRFKQLLNPAQ
 VFDLMNGGPHLGLRFLQKFDTRILVCGGDGSGVWVLEIDSLNLHKQCQLGVLPLGTGNDLARVLGWGS
 ACDDDTQLPQILEKLERASTKMLDRWSVMAYEAKLPRQASSSTVTEDFSEDSEVQQLFYEDSVA AHL SK
 ILTSDQHSSVVISSAKVLCETVKDFVARVKGAYEKTTESESESEVMAKKCSVLKEKLD SLLKTLDDDESQAS
 SSSLNPPPTIAEEAEDGDGSGSICGSTGDRLVASACPARPQIFRPREQLMLRANSLKKAIRQIIIEHTEKA
 VDEQNAQTQE QEGFVLGLSEEEKMDHRVCPPLSHSESGVPGKGRSQRKVS KSPCEKLISKGSLSGSSA
 SLPPQPGSRDGLPALNTKILYPNVRAGMSGSLPGG SVISRLLINADPFNSEPETLEYYTEKVMN NYFGI
 GLDAKISLDFNKRDEHPEKCRSRTKNMMWYGVLTGKELLHRTYKNLEQKVLLECDGRPIPLPSLQGI AV
 LNIPSYAGGTNFWGGTKEDDTFAAPSFDDKILEVAVFGSMQMAVSRVIRLQHHRIAQCRTVKISILGDE
 GVPVQVDGEAWVQPPGYIRIVHKNRAQTLTRDRAFESTLKS WEDKQKCELP RPPSCSLHPEMLSEEEATQ
 MDQFGQAAGVLIHSIREIAQSHRDMEQELAHAVNASSKSM DRVYGKPR TTEGLNCSFVLEMVNNFRALRS
 ETELLLSGKMALQLDPPQKEQLGSALAEMDRQLRRLADTPWLCQSAEPGDEESVMLDLAKRSRSGKFR LV
 TKFKKEKNNKNKEAHSSLGAPVHLWGTEEVAAWLEHLSLCEYKDIFTRHDIRGSELLHLERRDLKDLGVT
 KVGHMKRILCGIKELSRSA PAVEA


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK
Predicted MW: 134.3 kDa
Concentration: >0.05 µg/µL as determined by microplate BCA method



[View online »](#)

Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_690618
Locus ID:	8527
UniProt ID:	Q16760
RefSeq Size:	6294
Cytogenetics:	2q37.1
RefSeq ORF:	3642
Synonyms:	DGK-delta; dgkd-2; DGKdelta
Summary:	This gene encodes a cytoplasmic enzyme that phosphorylates diacylglycerol to produce phosphatidic acid. Diacylglycerol and phosphatidic acid are two lipids that act as second messengers in signaling cascades. Their cellular concentrations are regulated by the encoded protein, and so it is thought to play an important role in cellular signal transduction. Alternative splicing results in two transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]
Protein Families:	Druggable Genome
Protein Pathways:	Glycerolipid metabolism, Glycerophospholipid metabolism, Metabolic pathways, Phosphatidylinositol signaling system

Product images:

Coomassie blue staining of purified DGKD protein (Cat# [TP317053]). The protein was produced from HEK293T cells transfected with DGKD cDNA clone (Cat# [RC217053]) using MegaTran 2.0 (Cat# [TT210002]).