

Product datasheet for TP317098

FGFR2 (NM_000141) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human fibroblast growth factor receptor 2 (FGFR2), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC217098 representing NM_000141 Red =Cloning site Green =Tags(s)

MVSWGRFICLVVWTMATLSLARPSFSLVEDTTLEPEEPPTKYQISQPEVYVAAPGESLEVRCLLKDAAVI
SWTKDGVHLGPNNRVTLIGEYLQIKGATPRDSGLYACTASRTVDSETWYFMVNVTDAISSGDEDDTDGA
EDFVSENSNNKRAPYWTNTEKMEKRLHAVPAANTVKFRCPAGGNPMPTMRWLKNGKEFKQEHRIGGY
KVR
NQHWSLIMESVWPSDKGNYTCVENEYGSINHTYHLDVVERSHPHPILQAGLPANASTVVGDDVEFVCKV
YSDAQPHIQWIKHVKEKNGSKYGPDGLPYLKVLAAGVNTTDKEIEVLYIRNVTFEDAGEYTCLAGNSIGI
SFHSAWLTVLPAPGREKEITASPDYLEIAIYICIGVFLIACMVVTVILCRMKNTTKKPDFSSQPAVHKLTK
RIPLRRQVTVSAESSMNSNTPLVRITRRLSSTADTPMLAGVSEYELPEDPKWEFPRDKLTLGKPLGEG
CFGQVMAEAVGIDKDKPKEAVTVAVKMLKDDATEKDLSDLVSEMEMMKMIGKHKNIINLLGACTQDGP
L
YVIVEYASKGNLREYLRARRPPGMEYSYDINRVPEEQMTFKDLVSCTYQLARGMEYLASQKCIHRDLAAR
NVLVTENNVMKIADFLGARDINNIDYKKTNGRLPVKWMapeALFDRVYTHQSDVWSFGVLMWEIFTL
G
GSPYPGIPVEELFKLLKEGHRMDKPANCTNELYMMMRDCWHAVPSQRPTFKQLVEDLDRILTLTTNEEYL
DLSQPLEQYSPSPDTRSSCSSGDDSVFSPDPMPYEPCLPQYPHINGSVKT


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	89.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol



[View online »](#)

Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_000132
Locus ID:	2263
UniProt ID:	P21802
RefSeq Size:	4654
Cytogenetics:	10q26.13
RefSeq ORF:	2463
Synonyms:	BBDS; BEK; BFR-1; CD332; CEK3; CFD1; ECT1; JWS; K-SAM; KGFR; TK14; TK25
Summary:	<p>The protein encoded by this gene is a member of the fibroblast growth factor receptor family, where amino acid sequence is highly conserved between members and throughout evolution. FGFR family members differ from one another in their ligand affinities and tissue distribution. A full-length representative protein consists of an extracellular region, composed of three immunoglobulin-like domains, a single hydrophobic membrane-spanning segment and a cytoplasmic tyrosine kinase domain. The extracellular portion of the protein interacts with fibroblast growth factors, setting in motion a cascade of downstream signals, ultimately influencing mitogenesis and differentiation. This particular family member is a high-affinity receptor for acidic, basic and/or keratinocyte growth factor, depending on the isoform. Mutations in this gene are associated with Crouzon syndrome, Pfeiffer syndrome, Craniosynostosis, Apert syndrome, Jackson-Weiss syndrome, Beare-Stevenson cutis gyrata syndrome, Saethre-Chotzen syndrome, and syndromic craniosynostosis. Multiple alternatively spliced transcript variants encoding different isoforms have been noted for this gene. [provided by RefSeq, Jan 2009]</p>
Protein Families:	Druggable Genome, Protein Kinase, Secreted Protein, Transmembrane
Protein Pathways:	Endocytosis, MAPK signaling pathway, Pathways in cancer, Prostate cancer, Regulation of actin cytoskeleton

Product images:

Coomassie blue staining of purified FGFR2 protein (Cat# TP317098). The protein was produced from HEK293T cells transfected with FGFR2 cDNA clone (Cat# [RC217098]) using MegaTran 2.0 (Cat# [TT210002]).