

Product datasheet for TP317222


TMEM190 (NM_139172) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human transmembrane protein 190 (TMEM190), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC217222 protein sequence Red =Cloning site Green =Tags(s) MLGCGIPALGLLLLLQGSADGNGIQGFFYPWSCEGDIWDRESCGGQAAIDSPNLCLRLRCCYRNGVCYHQ RPDENVRRKHMWALVWTCGLLLLSCSICLFWWAKRRDVLHMPGFLAGPCDMSKSVSLLSKHRGTTKT PS TGSVPVALSKESRDVEGGTEGEGTEEGEETEGEEEEED TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	19.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_631911</u>
Locus ID:	147744
UniProt ID:	<u>Q8WZ59</u>


[View online »](#)

RefSeq Size:	593
Cytogenetics:	19q13.42
RefSeq ORF:	531
Synonyms:	MDAC1
Protein Families:	Transmembrane

Product images:


Coomassie blue staining of purified TMEM190 protein (Cat# TP317222). The protein was produced from HEK293T cells transfected with TMEM190 cDNA clone (Cat# [RC217222]) using MegaTran 2.0 (Cat# [TT210002]).