

Product datasheet for TP317285

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

BAIAP2 (NM 017451) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Purified recombinant protein of Homo sapiens BAl1-associated protein 2 (BAIAP2), transcript

variant 2, 20 µg

Species: Human
Expression Host: HEK293T

Expression cDNA Clone >RC217285 representing NM_017451 or AA Sequence: Red=Cloning site Green=Tags(s)

Red=Cloning site Green=Tags(s)

MSLSRSEEMHRLTENVYKTIMEQFNPSLRNFIAMGKNYEKALAGVTYAAKGYFDALVKMGELASESQGSK ELGDVLFQMAEVHRQIQNQLEEMLKSFHNELLTQLEQKVELDSRYLSAALKKYQTEQRSKGDALDKCQAE LKKLRKKSQGSKNPQKYSDKELQYIDAISNKQGELENYVSDGYKTALTEERRRFCFLVEKQCAVAKNSAA YHSKGKELLAQKLPLWQQACADPSKIPERAVQLMQQVASNGATLPSALSASKSNLVISDPIPGAKPLPVP PELAPFVGRMSAQESTPIMNGVTGPDGEDYSPWADRKAAQPKSLSPPQSQSKLSDSYSNTLPVRKSVTPK NSYATTENKTLPRSSSMAAGLERNGRMRVKAIFSHAAGDNSTLLSFKEGDLITLLVPEARDGWHYGESEK TKMRGWFPFSYTRVLDSDGSDRLHMSLQQGKSSSTGNLLDKDDLAIPPPDYGAASRAFPAQTASGFKQRP

YSVAVPAFSQGLDDYGARSMSRNPFAHVQLKPTVTNDRCDLSAQGPEGREHGDGSARTLAGR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK
Predicted MW: 60.7 kDa

Concentration: $>0.05 \mu g/\mu L$ as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional

chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.

Storage: Store at -80°C.





BAIAP2 (NM_017451) Human Recombinant Protein - TP317285

Stability: Stable for 12 months from the date of receipt of the product under proper storage and

handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: NP 059345

Locus ID: 10458 UniProt ID: Q9UQB8

RefSeq Size: 2877

Cytogenetics: 17q25.3 RefSeq ORF: 1656

Synonyms: BAP2; FLAF3; IRSP53; WAML

Summary: The protein encoded by this gene has been identified as a brain-specific angiogenesis inhibitor


(BAI1)-binding protein. This adaptor protein links membrane bound G-proteins to cytoplasmic effector proteins. This protein functions as an insulin receptor tyrosine kinase substrate and suggests a role for insulin in the central nervous system. It also associates with a downstream effector of Rho small G proteins, which is associated with the formation of stress fibers and cytokinesis. This protein is involved in lamellipodia and filopodia formation in motile cells and may affect neuronal growth-cone guidance. This protein has also been identified as interacting with the dentatorubral-pallidoluysian atrophy gene, which is associated with an autosomal dominant neurodegenerative disease. Alternative splicing results in multiple transcript

variants encoding distinct isoforms.[provided by RefSeq, Jan 2009]

Protein Families: Druggable Genome

Protein Pathways: Adherens junction, Regulation of actin cytoskeleton

Product images:



Coomassie blue staining of purified BAIAP2 protein (Cat# TP317285). The protein was produced from HEK293T cells transfected with BAIAP2 cDNA clone (Cat# [RC217285]) using MegaTran 2.0 (Cat# [TT210002]).