

Product datasheet for TP317289M

TIM 1 (HAVCR1) (NM_001099414) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human hepatitis A virus cellular receptor 1 (HAVCR1), transcript variant 2, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC217289 representing NM_001099414 Red =Cloning site Green =Tags(s)
	MHPQWILSLILHLADSVAGSVKVGGEAGPSVTLPCHYSGAVTSMCWNRGSCSLFTCQNGIWTNGTHVT YRKDTRYKLLGDLSRRDVSILTIENTAVSDSGVYCCRVEHRGWFNDMKITVSLEIVPPKVTTPVITVPT VTVRTSTTVPTTTTVPMTTVPTTTVPTTMSIPTTTTTLTTMTVSTTTSVPTTTSIPTTTSVPVTTTST FVPPMPLPRQNHEPVATSPSSPQPAETHPTTLQGAIRREPTSSPLSYTTDGNDDTVESSDGLWNNNQDQ LFLEHSLLTANTTKGIYAGVCISVLVLLALLGVIIAKKYFFKKEVQQLSVSFSSLQIKALQNAVEKEVQA EDNIYIENSLYATD
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	37.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001092884</u>




[View online »](#)

Locus ID:	26762
UniProt ID:	Q96D42
RefSeq Size:	1493
Cytogenetics:	5q33.3
RefSeq ORF:	1092
Synonyms:	HAVCR; HAVCR-1; KIM-1; KIM1; TIM; TIM-1; TIM1; TIMD-1; TIMD1

Summary: The protein encoded by this gene is a membrane receptor for both human hepatitis A virus (HHAV) and TIMD4. The encoded protein may be involved in the moderation of asthma and allergic diseases. The reference genome represents an allele that retains a MTTVP amino acid segment that confers protection against atopy in HHAV seropositive individuals. Alternative splicing of this gene results in multiple transcript variants. Related pseudogenes have been identified on chromosomes 4, 12 and 19. [provided by RefSeq, Apr 2015]

Protein Families: Druggable Genome, Transmembrane

Product images:



Coomassie blue staining of purified HAVCR1 protein (Cat# [TP317289]). The protein was produced from HEK293T cells transfected with HAVCR1 cDNA clone (Cat# [RC217289]) using MegaTran 2.0 (Cat# [TT210002]).