

Product datasheet for TP317326M

BARHL2 (NM_020063) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human BarH-like homeobox 2 (BARHL2), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC217326 protein sequence Red =Cloning site Green =Tags(s) MTMEGASGSSFGIDTILSSASSGSPGMMNGDFRPLGEARTADFRSQATPSPCSEIDTVGTAPSSPISVTM EPPEHLVADATQHHLHLHHSQQPPPPAAAPTQSLQPLPQQQQLPPQQPPPPPPQQLGSAASAPRT STS SFLIKDILGDSKPLAACAPYSTSVSSPHHTPKQESNAVHESFRPKLEQEDSKTKLDKREDSQSDIKCHGT KEEGDREITSSRESPPVRAKKPRKARTAFSDHQLNQLERSFERQKYSVQDRMDLAAALNLDTQVKTWY QNRRTKWKRQTAVGLELLAEAGNYSALQRMFSPYFYHPSLLGSMDSSTAAAAAAAMYSSMYRTPPAHP P QLQRPLVPRVLHGLGPGGQPALNPLSSPIPGTPHPR TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	41.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.


[View online »](#)

RefSeq:	<u>NP_064447</u>
Locus ID:	343472
UniProt ID:	<u>Q9NY43</u>
RefSeq Size:	1496
Cytogenetics:	1p22.2
RefSeq ORF:	1161
Summary:	Potential regulator of neural basic helix-loop-helix genes.[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified BARHL2 protein (Cat# [TP317326]). The protein was produced from HEK293T cells transfected with BARHL2 cDNA clone (Cat# [RC217326]) using MegaTran 2.0 (Cat# [TT210002]).