

Product datasheet for **TP317359**

PSMD10 (NM_170750) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human proteasome (prosome, macropain) 26S subunit, non-ATPase, 10 (PSMD10), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC217359 representing NM_170750 Red =Cloning site Green =Tags(s)
	 MEGCVSNLMVCNLAYSGKLEELKESILADKSLATRTDQDSRTALHWACCSAGHTEIVEFLLQLGVPVNDKD DAGWSPLHIAASAGRDEIVKALLGKGAQVNAVQNQNGCTPLHYAASKNRHEIAVMLLEGGANPDAKDHYEA TAMHRAAAKDT TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	16 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_736606</u>
Locus ID:	5716
UniProt ID:	<u>O75832</u>



[View online »](#)

RefSeq Size: 1467


Cytogenetics: Xq22.3

RefSeq ORF: 453

Synonyms: dj889N15.2; p28; p28(GANK)

Summary: This gene encodes a subunit of the PA700/19S complex, which is the regulatory component of the 26S proteasome. The 26S proteasome complex is required for ubiquitin-dependent protein degradation. This protein is a non-ATPase subunit that may be involved in protein-protein interactions. Aberrant expression of this gene may play a role in tumorigenesis. Two transcripts encoding different isoforms have been described. Pseudogenes have been identified on chromosomes 3 and 20.[provided by RefSeq, Mar 2011]

Product images:



Coomassie blue staining of purified PSMD10 protein (Cat# TP317359). The protein was produced from HEK293T cells transfected with PSMD10 cDNA clone (Cat# [RC217359]) using MegaTran 2.0 (Cat# [TT210002]).