

Product datasheet for **TP317380M**

FAM156B (NM_001099684) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens family with sequence similarity 156, member B (FAM156B), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC217380 representing NM_001099684 Red =Cloning site Green =Tags(s)
	MDPLQKRNPASPSKSSPMTAAETSQEGPAPSQPSYSEQPMMGLSNLSPGPGHSQAVPLPEGLLHQRYR EE KTLEERRWERLEFLQRKKAFLRHVRRRHRDHMAPYAVGREARISPLGDRSQNRQCECRYCQSHRPNLS G IPGESNRAPHPSWETLVQGLSGLTSLGTNQPGLPEAALQPQETEEKRQRERQQESKIMFQRLKQWL EEN TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	24.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001093154</u>



[View online »](#)

Locus ID: 727866
 UniProt ID: [Q8NDB6](#)
 RefSeq Size: 1701
 Cytogenetics: Xp11.22
 RefSeq ORF: 639
 Synonyms: TMEM29B

Product images:



Coomassie blue staining of purified FAM156B protein (Cat# [TP317380]). The protein was produced from HEK293T cells transfected with FAM156B cDNA clone (Cat# [RC217380]) using MegaTran 2.0 (Cat# [TT210002]).