

Product datasheet for TP317538

IFRD1 (NM_001550) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human interferon-related developmental regulator 1 (IFRD1), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC217538 protein sequence Red =Cloning site Green =Tags(s)

MPKNKKRNTPHRGSSAGGGGSGAAAATAATAGGQHRNVQPFSDDEDASIETMSHCSGYSDPSSFAEDGPEV
LDEEGTQEDLEYKLGKGLIDLTDKSAKTRQAALEGIKNALASKMLYEFILERRMTLTDSIERCLKKKGKSD
EQRAAAALASVLCIQLGPGIESEEILKTLGPILKKIICDGSASMQRQTCATCFGVCCFIATDDITELYS
TLECLENIFTKSYLKEKDTTVICSTPNTVLHISLLAWTLLLTICPINEVKKKLEMHFHKLPSLLSCDDV
NMRIAAGESLALLFELARGIESDFFYEDMESLTQMLRALATDGNKHRKVDRKQKRSVFRDVLRAVEERD
FPTETIKFGPERMYIDCWVKKHTYDTFKEVLGSGMQYHLQSNEFLRNVFELGPPVMLDAATLKTMKISRF
ERHLYNSAAFKARTKARSKCRDKRADVGEFF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	50.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP_001541](#)

Locus ID: 3475

UniProt ID: [Q00458](#), [A4D0U1](#)

RefSeq Size: 3305


Cytogenetics: 7q31.1

RefSeq ORF: 1353

Synonyms: PC4; TIS7

Summary: This gene is an immediate early gene that encodes a protein related to interferon-gamma. This protein may function as a transcriptional co-activator/repressor that controls the growth and differentiation of specific cell types during embryonic development and tissue regeneration. Mutations in this gene are associated with sensory/motor neuropathy with ataxia. This gene may also be involved in modulating the pathogenesis of cystic fibrosis lung disease. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Oct 2010]

Product images:



Coomassie blue staining of purified IFRD1 protein (Cat# TP317538). The protein was produced from HEK293T cells transfected with IFRD1 cDNA clone (Cat# [RC217538]) using MegaTran 2.0 (Cat# [TT210002]).