

## Product datasheet for **TP317663**

### SCGB1C1 (NM\_145651) Human Recombinant Protein

#### Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human secretoglobin, family 1C, member 1 (SCGB1C1), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC217663 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MKGSRALLLVALTLFCICRMATGEDNDEFFMDFLQTLVGTPEELYEGTLGKYNVNEDAKAAMTELKSCR DGLQPMHKAELVKLLVQVLGSQDGA
	<b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	10.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_663626</a>
Locus ID:	147199
UniProt ID:	<a href="#">Q8TD33</a>
RefSeq Size:	431
Cytogenetics:	11p15.5



[View online »](#)

RefSeq ORF: 285  
Synonyms: RYD5  
Protein Families: Secreted Protein

### Product images:



Coomassie blue staining of purified SCGB1C1 protein (Cat# TP317663). The protein was produced from HEK293T cells transfected with SCGB1C1 cDNA clone (Cat# [RC217663]) using MegaTran 2.0 (Cat# [TT210002]).