

## Product datasheet for **TP317696M**

### SIAT4A (ST3GAL1) (NM\_003033) Human Recombinant Protein

#### Product data:

|                                       |   |
|---------------------------------------|---|
| Product Type:                         | Recombinant Proteins  |
| Description:                          | Recombinant protein of human ST3 beta-galactoside alpha-2,3-sialyltransferase 1 (ST3GAL1), transcript variant 1, 100 µg |
| Species:                              | Human   |
| Expression Host:                      | HEK293T   |
| Expression cDNA Clone or AA Sequence: | >RC217696 protein sequence<br><b>Red</b> =Cloning site <b>Green</b> =Tags(s)  |

MVTLRKRTLKVLTFVLVLFILTSFFLNYSHTMVATTWFPKQMVLELSENLRKLIKHRPCTCTHCIGQRKL  
SAWFDERFNQTMQPLLTAQNALLEDDTYRWWLRLQREKKPNNLNDTIKELFRVVPGNVDPMLEKRSVGCR  
RCAVVGNSGNLRESSYGPEIDSHDFVLRMNKAPTAGFEADVGTKTTHHLVYPESFRELGDNVSMILVVPFK  
TIDLEWVVSAITTGTISHTYIPVPAKIRVKQDKILYHPAFIKYVFDNWLQGHGRYPSTGILSVIFSMHV  
CDEVDLYGFGADSKGNWHHYWENNPAGAFRKTGVHDADFESNVTATLASINKIRIFKGR

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**


|                |  |
|----------------|--|
| Tag:           | C-Myc/DDK  |
| Predicted MW:  | 38.9 kDa   |
| Concentration: | >0.05 µg/µL as determined by microplate BCA method   |
| Purity:        | > 80% as determined by SDS-PAGE and Coomassie blue staining  |
| Buffer:        | 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol   |
| Preparation:   | Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.                                     |
| Note:          | For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process. |
| Storage:       | Store at -80°C.  |
| Stability:     | Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.        |
| RefSeq:        | <a href="#">NP_003024</a>  |
| Locus ID:      | 6482   |



[View online »](#)

|                   |   |
|-------------------|---|
| UniProt ID:       | <a href="#">Q11201</a> , <a href="#">A0A024R9L6</a>   |
| RefSeq Size:      | 6971  |
| Cytogenetics:     | 8q24.22   |
| RefSeq ORF:       | 1020  |
| Synonyms:         | 1; Gal-NAc6S; SIAT4A; SIATFL; ST3GalA; ST3GalA.1; ST3GalIA; ST3O  |
| Summary:          | <p>The protein encoded by this gene is a type II membrane protein that catalyzes the transfer of sialic acid from CMP-sialic acid to galactose-containing substrates. The encoded protein is normally found in the Golgi but can be proteolytically processed to a soluble form. Correct glycosylation of the encoded protein may be critical to its sialyltransferase activity. This protein, which is a member of glycosyltransferase family 29, can use the same acceptor substrates as does sialyltransferase 4B. Two transcript variants encoding the same protein have been found for this gene. Other transcript variants may exist, but have not been fully characterized yet. [provided by RefSeq, Jul 2008]</p> |
| Protein Families: | Secreted Protein, Transmembrane   |
| Protein Pathways: | Glycosphingolipid biosynthesis - ganglio series, Glycosphingolipid biosynthesis - globo series, Keratan sulfate biosynthesis, Metabolic pathways, O-Glycan biosynthesis   |

### Product images:



Coomassie blue staining of purified ST3GAL1 protein (Cat# [TP317696]). The protein was produced from HEK293T cells transfected with ST3GAL1 cDNA clone (Cat# [RC217696]) using MegaTran 2.0 (Cat# [TT210002]).